



A


Simax® Büchner Flasks

SIMAX®

- Heavy wall flask
- With side hose connection for vacuum filtration

- Borosilicate 3.3 glass

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA38000	2420/250	250	145 x 85	35	10	£283.00
FLA38020	2420/500	500	175 x 105	35	8	£268.00
FLA38040	2420/1000	1000	230 x 135	45	8	£414.00
FLA3806	2420/2000	2000	255 x 165	60	Each	£54.10



A

PYREX® Büchner Flasks

PYREX®

- Büchner type flasks for vacuum filtration
- Heavy walls provide increased mechanical strength
- 'Studded' base to minimise surface abrasion
- Integral glass side-arm to accept 9mm bore flexible tubing
- Complies with ISO 6556

Code	Alt Ref	Capacity, mL	Height x Ø, mm	Neck I.D., mm	Pack	Price
FLA37000	1170/02M	100	105 x 64	24	10	£226.00
FLA37020	1170/04M	250	155 x 85	35	10	£263.00
FLA37040	1170/06M	500	185 x 105	35	10	£294.00
FLA37060	1170/08M	1000	230 x 135	45	5	£188.00
FLA3708	1170/10D	2000	260 x 166	60	Each	£99.00
FLA3707	1170/12M	3000	335 x 187	68	5	£715.00
FLA3710	1170/14D	5000	385 x 220	68	Each	£189.00

GO SAFE TRUST IN DWK

Safety innovations
and expertise for your
daily lab work



Explore more at
scientificlabs.co.uk

DWK
LIFE SCIENCES