

DURAN® Büchner Flasks



- DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high performance under vacuum giving long service
- Very resistant to chemical attack
- Suitable for autoclave use
- With marking label
- Available with glass or plastic side arm
- Capacities up to 2L are Erlenmeyer (conical) form to ISO 6556:2012; 3L and above are bottle form

Glass Sidearm

Code	Alt Ref	Capacity, mL	Height, mm	Neck Ø, mm	Pack	Price
FLA39000	212012409	100	105	24	10	£247.00
FLA39020	212013602	250	155	34	10	£288.00
FLA39040	212014401	500	185	34	10	£320.00
FLA39060	212015406	1000	230	45	10	£457.00
FLA3908	212016308	2000	255	60	Each	£110.00
FLA3910	211916802	3000	295	58	Each	£119.00
FLA3912	211917301	5000	360	68	Each	£200.00
FLA3914	211918606	10000	420	70	Each	£295.00
FLA3916	211918803	15000	500	70	Each	£375.00
FLA3918	211919105	20000	535	70	Each	£468.00



Plastic Sidearm

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA39200	212042452	100	105 x 64	24	10	£197.00
FLA39220	212043654	250	420 x 240	34	10	£254.00
FLA39240	212044453	500	155 x 85	34	10	£320.00
FLA39260	212045458	1000	500 x 257	45	10	£477.00
FLA3928	212046351	2000	185 x 105	60	Each	£99.00
FLA3927	211946854	3000	535 x 290	58	Each	£123.00
FLA3929	211947353	5000	230 x 135	68	Each	£171.00
FLA3921	211948658	10000	295 x 170	70	Each	£271.00
FLA3925	211949157	20000	255 x 166	70	Each	£441.00
FLA3923	211948855	15000	360 x 185	70	Each	£354.00



Plastic Sidearm Spares

Code	Alt Ref	Description	Pack	Price
FLA3940	292585407	KECK™ hose connection assembly set	10	£56.00
FLA3944	292201009	Spare silicon seal	10	£12.60