

Sartorius Secura® Precision Balances

SARTORIUS

Reduce your risk with integrated and intelligent protection systems. Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, built in protection systems are included for complete reliability and regulatory compliance – the safe and secure way. Real time guidance prompts for levelling, automatic internal adjustment, and 100% traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Integrated protection systems for minimising your risk of obtaining incorrect results
- Real time alarm messaging and interactive user guidance for levelling
- Automatic detection and active monitoring to ensure you are within the acceptable USP operating range for reliable results
- Fully automatic internal adjustment function, isoCAL, for the highest possible accuracy at all times
- Cal Audit Trail for easy, traceable documentation of all calibration and adjustment data
- Safe and clear documentation by sample and batch identifiers for keeping perfect track of weighing data
- GLP/GMP compliant printouts
- Standard password security features for protection against tampering and unauthorised access




Code	Alt Ref	Capacity, g	Readability, g	Pan Ø, mm	Windshield	Dims, w x d x h, mm	Price
BAL7008	SECURA213-1S	210	0.001	120	Manual	216 x 360 x 320	£3034.00
BAL7012	SECURA313-1S	310	0.001	120	Manual	216 x 360 x 320	£3124.00
BAL7016	SECURA513-1S	510	0.001	120	Manual	216 x 360 x 320	£3912.00
BAL7020	SECURA1102-1S	1100	0.01	180	No	216 x 360 x 95	£2723.00
BAL7024	SECURA2102-1S	2100	0.01	180	No	216 x 360 x 95	£3000.00
BAL7028	SECURA3102-1S	3100	0.01	180	No	216 x 360 x 95	£3156.00
BAL7032	SECURA5102-1S	5100	0.01	180	No	216 x 360 x 95	£4067.00

Sartorius Quintix® Precision Balances

SARTORIUS

Quintix® makes laboratory work more enjoyable. The intuitive touchscreen in combination with the most precise weighing technology make Quintix® laboratory balances the only choice to achieve the best results. Benefit from improved workflows by easy to use application programs. If you need to generate GLP compliant documents, you will save considerable time and effort with "plug and work" connectivity.

- Intuitive operation of the balance without having to read a manual thanks to the new user interface and colour touch screen
- Integrated application programs and intuitive user guidance for greater efficiency in your daily lab routine
- Highest accuracy and convenience ensured at all times by the isoCAL fully automatic internal adjustment feature



Code	Alt Ref	Capacity, g	Readability, g	Pan Ø, mm	Windshield	Dims, w x d x h, mm	Price
BAL7044	QUINTIX213-1S	210	0.001	120	Manual	216 x 360 x 320	£2512.00
BAL7048	QUINTIX313-1S	310	0.001	120	Manual	216 x 360 x 320	£2701.00
BAL7052	QUINTIX513-1S	510	0.001	120	Manual	216 x 360 x 320	£3223.00
BAL7056	QUINTIX612-1S	610	0.01	180	No	216 x 360 x 95	£2001.00
BAL7060	QUINTIX1102-1S	1100	0.01	180	No	216 x 360 x 95	£2256.00
BAL7064	QUINTIX2102-1S	2100	0.01	180	No	216 x 360 x 95	£2549.00
BAL7068	QUINTIX3102-1S	3100	0.01	180	No	216 x 360 x 95	£2645.00
BAL7080	QUINTIX5102-1S	5100	0.01	180	No	216 x 360 x 95	£3400.00
BAL7076	QUINTIX5101-1S	5100	0.1	180	No	216 x 360 x 95	£2301.00
BAL7072	QUINTIX5100-1S	5100	1.0	180	No	216 x 360 x 95	£1868.00