



A

DURAN® Wide Mouth Erlenmeyer Flasks



- ISO 24450:2005
- DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass Class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service
- Very resistant to chemical attack
- Suitable for hotplate and autoclave use
- Marking label, graduations and retrace code which allows downloading of individual batch certificate

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Pack	Price
FLA42200	212261405	25	70 x 43	31	10	£55.00
FLA42220	212261705	50	85 x 51	34	10	£55.00
FLA42240	212262401	100	105 x 64	34	10	£56.00
FLA42260	212263209	200	131 x 79	50	10	£58.00
FLA42280	212263603	250	140 x 85	50	10	£62.00
FLA42300	212263903	300	156 x 87	50	10	£71.00
FLA42320	212264402	500	175 x 105	50	10	£80.00
FLA42340	212265407	1000	220 x 131	50	10	£128.00
FLA4236	212266309	2000	276 x 153	72	Each	£23.90
FLA4238	212276807	3000	285 x 190	106	Each	£81.00
FLA4239	212277306	5000*	322 x 220	108	Each	£105.00
FLA4235	212278602	10000*	420 x 285	147	Each	£208.00

* Flasks with a capacity >3L are without graduations or retrace code



A

DURAN® Narrow Neck Erlenmeyer Flasks, with Screw Caps



- ISO 1773:1997
- DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass Class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service
- Very resistant to chemical attack
- With marking label, graduations and retrace code which allows downloading of individual batch certificate
- Supplied with PBT caps
- Also compatible with DURAN® autoclaving caps BOT2103/BOT2105

Supplied in packs of 10

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Neck Ø, mm	Thread	Price
FLA40420	218032451	100	109 x 64	25	GL25	£309.00
FLA40440	218033653	250	149 x 85	32	GL32	£298.00
FLA40460	218034452	500	180 x 105	32	GL32	£314.00
FLA40480	218035457	1000	225 x 131	32	GL32	£472.00