


DURAN® Culture Flasks Erlenmeyer Rimless

- Straight rimless neck for metal caps
- DURAN® borosilicate glass 3.3 to ISO 3585:1998; glass Class 1 to USP, EP & JP
- Tightly controlled glass quality, wall thickness and geometry ensure high resistance to thermal shock giving long service
- Very resistant to chemical attack
- Suitable for hotplate and autoclave use
- With marking label, graduations and Retrace Code which allows downloading of individual batch certificate
- Neck Ø: 38mm
- Supplied without caps
- For DURAN® Erlenmeyer flasks with screw caps, see FLA40420 series

Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Pack	Price
FLA42700	217712401	100	114 x 64	10	£173.00
FLA42720	217713603	250	149 x 85	10	£192.00
FLA42740	217713903	300	161 x 87	10	£96.00
FLA42760	217714402	500	183 x 105	10	£119.00
FLA42780	217715407	1000	229 x 131	10	£247.00
FLA4280	217716309	2000	302 x 166	Each	£31.60




Rimless Flask Metal Caps

- Suitable for neck Ø 38mm
- Metal caps allow an exchange of air in the media bottle, creating a positive environment for the growth of microorganisms

Supplied in packs of 10

Code	Alt Ref	Material	Colour	Price
FLA4266	290122406	Stainless steel	Silver	£149.00
FLA4268	290132407	Aluminium	Blue	£100.00



Thermo Scientific™ Nalgene™ Erlenmeyer Flasks

Thermo
SCIENTIFIC

- Manufactured from polypropylene copolymer
- Excellent chemical compatibility and shatter-proof construction.
- Translucent wall makes liquid level easily visible
- 2000mL size has moulded-in graduations
- Mouth sized to accept standard sized stoppers (sold separately)
- Repeatedly autoclavable for long life and rigorous use

Code	Alt Ref	Capacity, mL	No. Stopper	Pack	Price
FLA4154	4102-0125	125	17	12	£102.00
FLA4156	4102-0250	250	25	6	£70.00
FLA4158	4102-0500	500	29	4	£67.00
FLA4160	4102-1000	1000	31	2	£52.70
FLA4162	4102-2000	2000	41	2	£64.00

