

PYREX® Volumetric Flasks

PYREX®

- Complies with ISO 1042 and DIN 12664
- Calibrated to the tolerance required by the United States Pharmacopoeia (USP) which is tighter than those of ISO/DIN standards
- Each flask is laser etched with an individual serial number or batch code and date of manufacture
- Works certified flasks are supplied with printed calibration certificate
- Conformity 'H' certified
- Manufactured from chemically resistant PYREX® glass
- Inscriptions in white enamel
- Polyethylene stopper



A

Class A

Code	Alt Ref	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA73300	3151/04	5	0.020	10/19	5	£45.20
FLA73320	3151/06	10	0.020	10/19	5	£45.20
FLA73340	3151/08	20	0.03	10/19	5	£53.10
FLA73360	3151/10	25	0.03	10/19	5	£53.20
FLA73380	3151/12	50	0.05	12/21	5	£55.00
FLA73400	3151/14	100	0.08	14/23	5	£68.00
FLA73420	3151/16	200	0.10	14/23	5	£92.00
FLA73440	3151/18	250	0.12	14/23	2	£37.50
FLA73460	3151/20	500	0.20	19/26	2	£50.50
FLA73480	3151/22	1000	0.30	24/29	2	£67.00
FLA7350	3151/24	2000	0.50	29/32	Each	£59.00
FLA7352	3151/26	5000	1.00	34/35	Each	£148.00

Class A, Works Certified

Code	Alt Ref	Capacity, mL	Tolerance, ± mL	Stopper Size	Pack	Price
FLA7101	3151/01WC	1	0.020	7/16	Each	£121.00
FLA7103	3151/02WC	2	0.020	7/16	Each	£121.00
FLA73000	3151/04WC	5	0.020	10/19	5	£79.00
FLA73020	3151/06WC	10	0.020	10/19	5	£79.00
FLA73040	3151/08WC	20	0.03	10/19	5	£87.00
FLA73060	3151/10WC	25	0.03	10/19	5	£86.00
FLA73080	3151/12WC	50	0.05	12/21	5	£88.00
FLA73100	3151/14WC	100	0.08	14/23	5	£100.00
FLA73120	3151/16WC	200	0.10	14/23	5	£129.00
FLA73140	3151/18WC	250	0.12	14/23	2	£49.40
FLA73160	3151/20WC	500	0.20	19/26	2	£64.00
FLA73180	3151/22WC	1000	0.30	24/29	2	£79.00
FLA7320	3151/24WC	2000	0.50	29/32	Each	£66.00
FLA7322	3151/26WC	5000	1.00	34/35	Each	£174.00