

Heathrow Scientific Lead Flask Rings



- For stabilising glassware and plasticware or for when immersed in water baths
- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- O-shaped

Code	Alt Ref	Fits Size, mL	Ext Ø, mm	Int Ø, mm	Price
FLA4900	HS8882A	125 - 500	77	48	£32.10
FLA4902	HS8882B	250 - 1000	90	51	£35.40
FLA4904	HS8882C	500 - 2000	98	57	£42.80
FLA4906	HS8882D	1000 - 2000	115	70	£55.00



Cork Flask Rings

- For supporting round bottom flasks
- Thickness: 45mm

Supplied in packs of 5

Code	Ext Ø, mm	Int Ø, mm	Price
COR1100	75	45	£33.50
COR1102	115	85	£51.23
COR1104	150	120	£66.00
COR1106	195	160	£74.27
COR1108	240	200	£110.00



Aluminium Foil

- Gauge 11µm ±1µm
- Does not harbour or promote bacterial growth

Code	Width, mm	Length, m	Price
ALU1005	300	75	£33.70
ALU1011	450	75	£43.61
ALU1015	500	90	£56.48



Aluminium Foil

- In dispenser box
- Dims, w x l: 326mm x 200m

Code	Thickness, µm	Price
39400	10	£21.00
39402	14	£20.10

