

## PYREX® Separating Funnels, Pear Shape

PYREX®

- Complies with ISO 4800 type 2
- Manufactured from PYREX® borosilicate glass for corrosion-free performance
- With interchangeable PTFE stopcock key
- Supplied with glass and polyethylene stopper

Code	Alt Ref	Capacity, mL	Stopper Size	Stopcock Bore, mm	Price
FUN4300	2180/02M	50	19/26	1.5	£74.00
FUN4322	2180/04M	100	19/26	1.5	£79.00
FUN4324	2180/06M	250	24/29	2.5	£82.00
FUN4306	2180/08M	500	24/29	4.0	£86.00
FUN4328	2180/10M	1000	29/32	4.0	£99.00
FUN4310	2180/12M	2000	29/32	6.0	£140.00
FUN4312	2180/14M	5000	34/35	6.0	£207.00



## PYREX® Glass Funnels

PYREX®

- Short stem addition funnels
- With 60° bowl angle for precise fitting of filter paper
- Manufactured from chemical and heat resistant PYREX® borosilicate glass
- Suitable for filtering a wide range of chemical solutions



Code	Alt Ref	Dims, h x rim Ø, mm	Stem Length, mm	Filter Paper Size, mm	Pack	Price
FUN11500	2140/01M	60 x 35	35	45 - 55	10	£116.00
FUN11520	2140/02M	78 x 45	45	55 - 70	10	£116.00
FUN11540	2140/03M	125 x 70	75	70 - 90	10	£126.00
FUN11560	2140/06M	144 x 80	82	90 - 150	10	£145.00
FUN11580	2140/10M	186 x 105	106	150 - 185	10	£151.00
FUN1160	2140/11M	220 x 125	122	185 - 240	Each	£19.80
FUN1162	2140/12M	270 x 150	154	240 - 270	Each	£30.30
FUN1164	2140/14M	345 x 200	193	320 - 385	Each	£67.00

## PYREX® Glass Powder Funnels

PYREX®

- Ideal for transferring powders
- With short wide stem for smooth flow
- Manufactured from chemically resistant PYREX® borosilicate glass

Code	Alt Ref	Top Ø, mm	Stem Length, mm	Pack	Price
FUN20000	2130/02M	80	25	10	£183.00

