



## Sartorius Practum® Precision Balances

SARTORIUS

Get the best value for your money, without compromises in precision and reliability. Rely on consistent readings and excellent repeatability. Practum® has everything you need for standard applications.


- Accurate, consistent readings with excellent repeatability based on Sartorius quality and technology designed in Germany
- High reliability and durability provided by mechanical overload protection
- Extra large feet and up front level indicator for easy levelling
- Supervisor lock for preventing unintentional changes to the setup settings
- Effortless weighing thanks to the ergonomic and sturdy draft shield (0.001g versions only)

Code	Alt Ref	Capacity, g	Readability, g	Pan Ø, mm	Windshield	Price
<b>BAL7088</b>	PRACTUM213-1S	210	0.001	120	Manual	£1120.00
<b>BAL7090</b>	PRACTUM313-1S	310	0.001	120	Manual	£1270.00
<b>BAL7092</b>	PRACTUM612-1S	610	0.01	180	No	£990.00
<b>BAL7094</b>	PRACTUM1102-1S	1100	0.01	180	No	£1070.00
<b>BAL7096</b>	PRACTUM2102-1S	2100	0.01	180	No	£1280.00
<b>BAL7100</b>	PRACTUM5101-1S	5100	0.1	180	No	£1070.00
<b>BAL7098</b>	PRACTUM5100-1S	5100	1.0	180	No	£809.00

## Sartorius Entris® II Essential Line Precision Balances

SARTORIUS

Every Sartorius balance offers quality, value and consistency. No matter what you're weighing, the new Entris® II is always the right choice. Offering unrivalled value at a budget price, this is the only balance in its class featuring isoCAL, LED touch technology and 12 built in applications. Backed by almost 150 years of German engineering expertise, and available in 40 different models, you'll easily find an Entris® II balance that exactly meets your specific weighing needs.



- Highly accurate results year in year out – guaranteed via Monolithic weigh cell technology invented by Sartorius
- Fastest stabilisation time in its class – using state of the art weighing sensors
- Optional isoCAL (Internal Calibration & Adjustment) for total assurance of accurate weighing results
- Hybrid screen for excellent readability and use, combines intuitive, wear resistant LED and touch technology
- 12 built in applications with GLP/GMP compliant printout/data output
- USB Type C and RS232 9 pin interfaces
- Password protection for secure operation

Code	Alt Ref	Capacity, g	Readability, mg	Pan Dims, mm	Windshield	Calibration	Price
<b>BAL1166</b>	BCE223-1S	220	1	Ø 120	Manual	External	£898.00
<b>BAL1168</b>	BCE323-1S	320	1	Ø 120	Manual	External	£990.00
<b>BAL1170</b>	BCE423-1S	420	1	Ø 120	Manual	External	£1090.00
<b>BAL1172</b>	BCE623-1S	620	1	Ø 120	Manual	External	£1280.00
<b>BAL1174</b>	BCE653-1S	650	1	Ø 120	Top access	External	£1470.00
<b>BAL1176</b>	BCE622-1S	620	10	182 x 182	—	External	£804.00