


Polypropylene Funnels

Code	Top Ø, mm	Pack	Price
FUN1250	40	10	£10.82
FUN1260	160	Each	£3.50



Hirsch Funnels

- Porcelain
- With perforated support plate

Code	Top Ø, mm	Plate Ø, mm	Price
FUN3020	50	18	£24.00
FUN3022	70	27	£18.00



PYREX® Büchner Funnels

PYREX®

- Büchner filter funnels with integral sintered glass discs
- Disc and funnel body manufactured from chemically resistant PYREX® borosilicate glass
- For convenient filtration of solutions
- No filter papers required
- Highly resistant to most corrosive reagents including ammonia and concentrated sulphuric acid
- Choice of diameters and porosities



Code	Alt Ref	Disc Ø, mm	Capacity, mL	Porosity	Height x Rim Ø, mm	Pack	Price
FUN71020	3620/02	30	35	1	160 x 38	5	£148.00
FUN71040	3620/04	30	35	2	160 x 38	5	£148.00
FUN71060	3620/06	30	35	3	160 x 38	5	£149.00
FUN71080	3620/08	30	35	4	160 x 38	5	£148.00
FUN71200	3630/01	40	50	0	170 x 54	5	£153.00
FUN71220	3630/02	40	50	1	170 x 54	5	£159.00
FUN71240	3630/04	40	50	2	170 x 54	5	£159.00
FUN71260	3630/06	40	50	3	170 x 54	5	£160.00
FUN71280	3630/08	40	50	4	170 x 54	5	£158.00
FUN71300	3630/10	40	50	5	170 x 54	5	£168.00
FUN71400	3700/01	60	125	0	180 x 71	5	£315.00
FUN71420	3700/02	60	125	1	180 x 71	5	£318.00
FUN71440	3700/04	60	125	2	180 x 71	5	£318.00
FUN71460	3700/06	60	125	3	180 x 71	5	£319.00
FUN71480	3700/08	60	125	4	180 x 71	5	£318.00
FUN71500	3700/10	60	125	5	180 x 71	5	£336.00
FUN71640	3710/04	95	500	2	250 x 106	5	£438.00
FUN71660	3710/06	95	500	3	250 x 106	5	£439.00
FUN71820	3720/02M	120	1000	1	290 x 135	2	£239.00
FUN71840	3720/04M	120	1000	2	290 x 135	2	£239.00
FUN71860	3720/06M	120	1000	3	290 x 135	2	£239.00
FUN71880	3720/08M	120	1000	4	290 x 135	2	£239.00
FUN7202	3722/02M	175	4000	1	435 x 190	Each	£264.00
FUN7204	3722/04M	175	4000	2	435 x 190	Each	£264.00
FUN7206	3722/06M	175	4000	3	435 x 190	Each	£269.00
FUN7208	3722/08M	175	4000	4	435 x 190	Each	£269.00