



Kango Polyurethane Stools



- Round stools in black polyurethane
- Chrome, 5-branch, scuff-resistant, reinforced base
- Chrome gas-lift column with height adjustment and shock absorber

Code	Alt Ref	Description	Glides/ Castors	Mechanism	Footring	Height, mm	Price
FUR5062	7NN10GHJR02905	Round stool	Nylon castors	Adjustable height	Yes	580 - 820	£284.00
FUR5362	7NN10GHJP02905	Round stool	Glides	Adjustable height	Yes	540 - 780	£206.00
FUR5063	7NN10GHJP00905	Round stool	Glides	Adjustable height	No	540 - 780	£202.00
FUR5328	7NS20GHJU00905	Saddle stool	Rubber castors	Adjustable height	No	580 - 820	£296.00
FUR5310	7NS 27GPJU 00 905	Saddle stool	Rubber castors	Adjustable height by foot control	No	600 - 700	£498.00



Kango Polyurethane Semi Stand Up Chairs



- Ergonomic seat in soft, black polyurethane, tilted forward
- Chrome, 5-branch, scuff-resistant, reinforced base with articulated glides
- Available with or without footrest
- Chrome gas-lift column with height adjustment from 550 to 790mm and shock absorber

Code	Alt Ref	Footring	Price
FUR5064	7NG20GHLPO0905	No	£289.00
FUR5372	7NG20GHLPO2905	Yes	£294.00

Labsit Chairs

Labsit is the new chair for all kinds of laboratories. Available with artificial leather or PU foam upholstery. Artificial leather is washable, low maintenance, and resistant to disinfectants while being soft and comfortable. It is PVC and phthalate free, resistant to blood and urine, has fire protection according to EN 1021 part 1 + 2, OEKO-TEX® standard 100 conform, and >50,000 Martindale rubs. PU foam upholstery is extremely robust, durable and washable. It is resistant to mechanical influences, disinfectants, blood and urine, antibacterial, PVC-free with fire protection according to EN 1021 part 1 + 2.

- Easy to clean with minimal joints, is resistant to disinfectants and solvents
- Anti-bacterial and anti-microbial properties
- Suitable for cleanrooms in air purity class 3, acc. to EN ISO 14644-1
- Low weight and integrated handle provide maximum mobility
- Compact structure to save space
- Comfortable and wear resistant
- Flex function replaces complex adjustments, the backrest and front edge of the seat follow the movements of the user and provides optimal support for the body in any position
- Meets the specifications of the GMP regulations
- Black polyamide base with glides and foot ring
- Seat dims, w x d: 450 x 420mm
- Backrest height: 420mm
- Ø 5 star base: 640mm



Labsit 2

- Laboratory chair with castors
- Seat height: 450 - 650mm

Code	Alt Ref	Upholstered Material	Seat Shell Colour	Price
FUR3500	9123-MG01-3277	Artificial leather	Blue	£217.00
FUR3502	9123-MG01-3285	Artificial leather	Anthracite	£188.00
FUR3504	9123-MG01-3279	Artificial leather	Orange	£188.00
FUR3506	9123-MG01-3280	Artificial leather	Green	£188.00
FUR3508	9123-2000-3277	Integral foam	Blue	£188.00
FUR3510	9123-2000-3285	Integral foam	Anthracite	£188.00
FUR3512	9123-2000-3279	Integral foam	Orange	£188.00
FUR3514	9123-2000-3280	Integral foam	Green	£188.00