

Tourniquets

- Designed for venipuncture or intravenous injection
- High quality
- Made from durable Velcro for easy use and fast release

Code	Pack	Price
HAE3902	10	£37.70



BD PAXgene™ Blood RNA Tubes



- Intended for immediate stabilisation of intracellular RNA, thereby yielding accurate and reproducible gene expression data
- Contains 6.9mL of proprietary RNA stabilisation additive
- Intracellular RNA is stabilised for up to 3 days at room temperature or 5 days at 2 to 8°C
- Can be frozen at temperatures of -20 to -70°C for long term storage
- Sterilisation: radiation, 10⁻⁶ SAL, ISO 11137-1
- Closure colour: red

Supplied in packs of 100



Code	Draw Volume, mL	Closure Type	Dims, h x Ø, mm	Pack	Price
762165	2.5	Hemogard™	100 x 16	100	£1382.00

BD Vacutainer® Serum Separator Tubes



- SST™ II *Advance* tubes for obtaining and separating a serum sample
- Contain Silica to activate clotting of the specimen, and a gel that forms a barrier between the clot and the serum after centrifugation
- Design enables centrifugation flexibility which provides increased laboratory workflow and process efficiency
- Mixing recommendations: gently invert to 180° and back 5 – 6 times
- Centrifuge instructions: 1300 – 2000g for 10 minutes at room temperature
- Closure colour: gold



Code	Capacity, mL	Dims, h x Ø, mm	Pack	Price
367956	3.5	75 x 13	1000	£205.00
367958	8.5	100 x 16	100	£35.00