

Epredia™ HM 430 Sliding Microtome Continued

epredia | Authorized Distributor

Accessories

Code	Alt Ref	Description	Price
HIS3034	705200	Disposable blade holder with eccentric fast clamping for low profile blades	£916.00
HIS3036	705340	Blade holder with eccentric fast clamping for high profile blades	£932.00
HIS3035	715480	Standard specimen clamp (45 x 80mm) for HM 430, HM 450	£607.00

Epredia™ HM 450 Electronic Sliding Microtome

epredia | Authorized Distributor

Includes automatic fine sectioning and trimming features with selectable specimen retraction. The integrated knife carrier design of the HM 450 Sliding Microtome ensures fast, comfortable and non-tiring operation.

- Stepping motor technology assures precise, accurate specimen advancement
- Trimming function allows for precision block trimming prior to the sectioning function
- User-friendly control panel allows for programming of settings and sectioning parameters
- Important sectioning information including section counter, step sum, remaining travel and specimen position memory function are retained
- Specimen fine orientation on X/Y-axes
- Knife carrier glides along a maintenance-free precision cross-roller bearing system
- Memory function for time-saving specimen advance to the knife
- Ergonomic handle for moving knife carrier
- Automatic or manual sectioning
- Note: requires standard specimen clamp (715480)



Specifications	
Section Thickness Range, µm	1 to 100
Vertical Feed Range, mm	40 max
Horizontal Blade Stroke, mm	190 max
Dims, w x d x h, mm	360 x 490 x 320

Code	Alt Ref	Supplied with	Price
HIS3014	910020	Steel knife holder, dust cover, operator manual	£8590.00

Epredia™ MX35 Ultra Disposable Microtome Blades

epredia | Authorized Distributor

Disposable blades featuring a proprietary steel hardening technology (Pink Technology) that delivers unbeatable longevity, making them excellent for cutting frozen sections.

- Unique angle profile that facilitates faster ribboning (34°)
- Low profile
- Proprietary process ensures a durable and consistent coating



Code	Alt Ref	Blade Angle, °	Pack	Price
HIS3050	3053835	34	50	£75.00