

Epredia™ HistoStar™ Embedding Workstations Continued

epredia | Authorized Distributor

Code	Alt Ref	Description	Price
HIS3350	A81000001	HistoStar embedding centre complete system	£13179.00
HIS3352	A81000004	HistoStar cold module	£5257.00
HIS3354	A81000007	HistoStar embedding module	£9766.00

Accessories

Code	Alt Ref	Description	Pack	Price
HIS3356	A81010030	Tamping tool kit, 1 x 20mm, 1 x 12mm	2	£40.10
HIS3358	A81010032	Magnifier assembly	Each	£285.00
HIS3360	A81010033	Heated forceps, red, 1mm tip	Each	£463.00
HIS3362	A81010034	Heated forceps, yellow, 2mm tip	Each	£463.00
HIS3364	A81010035	Heated forceps, blue, 4mm tip	Each	£465.00

Epredia™ CTM6 Coverslipper

epredia | Authorized Distributor

Get compact automation of routine glass coverslipping. With its compact footprint the CTM6 is the ideal instrument to free staff from the time-consuming task of routine coverslipping.

- Operator safety - loading and unloading of racks without opening the hood
- Internal vapour extraction system with fan and permanent air-flow avoids any leakage of vapours at openings
- Integrated charcoal filtration
- Bubble-free mountant refill, 100mL reservoir
- Tool-free cleaning and maintenance
- Easy conversion to other coverslip sizes - all required parts included



Code	Alt Ref	Dims, w x d x h, mm	Price
HIS3602	970010	275 x 410 x 350	£26943.00

Chemicals & Consumables



Laboratory consumables
and chemicals to
suit a wide range
of applications