

Simport® Scientific Histosette® II Cassettes



These patented cassettes keep specimens safely submerged in liquid and are resistant to the chemical action of most histological solvents. The efficient flow-through slots maximise fluid exchange and ensure proper drainage.

- The unobstructed anterior writing area is slanted at a 45° angle
- Made of acetal
- Suitable to be used with certain types of cassette labelling instruments
- Easy write-on marking surface for convenient identification



With Attached Lids

Lid is pre-attached to base. Back mounted locking device securely holds the lid in place, front hinge design allows single handed cassette manipulation.

Supplied in cases of 1500 (3 dispenser boxes of 500 cassettes)

Code	Alt Ref	Description	Colour	Price
HIS0065	M492-2	Tissue cassettes	White	£282.00
HIS0066	M492-3	Tissue cassettes	Pink	£282.00
HIS0067	M492-4	Tissue cassettes	Green	£276.00
HIS0068	M492-5	Tissue cassettes	Yellow	£282.00
HIS0069	M492-6	Tissue cassettes	Blue	£282.00
HIS0070	M492-7	Tissue cassettes	Peach	£282.00
HIS0071	M492-8	Tissue cassettes	Tan	£282.00
HIS0072	M492-9	Tissue cassettes	Grey	£282.00
HIS0073	M492-10	Tissue cassettes	Lilac	£282.00
HIS0074	M492-11	Tissue cassettes	Orange	£282.00
HIS0075	M492-12	Tissue cassettes	Aqua	£282.00

Code	Alt Ref	Description	Colour	Price
HIS0076	M493-2	Biopsy cassettes	White	£274.00
HIS0077	M493-3	Biopsy cassettes	Pink	£274.00
HIS0078	M493-4	Biopsy cassettes	Green	£274.00
HIS0079	M493-5	Biopsy cassettes	Yellow	£274.00
HIS0080	M493-6	Biopsy cassettes	Blue	£274.00
HIS0081	M493-7	Biopsy cassettes	Peach	£274.00
HIS0082	M493-8	Biopsy cassettes	Tan	£274.00
HIS0083	M493-9	Biopsy cassettes	Grey	£274.00
HIS0084	M493-10	Biopsy cassettes	Lilac	£274.00
HIS0085	M493-11	Biopsy cassettes	Orange	£274.00
HIS0086	M493-12	Biopsy cassettes	Aqua	£274.00



With Separate Lids

Cassettes and lids packaged separately.

Supplied in cases of 1000 (2 dispenser boxes of 500 cassettes, 1 dispenser box of 1000 lids)

Code	Alt Ref	Description	Colour	Price
HIS0087	M485-2	Tissue cassette	White	£142.00
HIS0088	M485-3	Tissue cassette	Pink	£142.00
HIS0089	M485-4	Tissue cassette	Green	£142.00
HIS0090	M485-5	Tissue cassette	Yellow	£142.00
HIS0091	M485-6	Tissue cassette	Blue	£142.00
HIS0092	M485-7	Tissue cassette	Peach	£142.00
HIS0093	M485-8	Tissue cassette	Tan	£142.00
HIS0094	M485-9	Tissue cassette	Grey	£142.00
HIS0095	M485-10	Tissue cassette	Lilac	£142.00
HIS0096	M485-11	Tissue cassette	Orange	£142.00
HIS0097	M485-12	Tissue cassette	Aqua	£142.00

Code	Alt Ref	Description	Colour	Price
HIS0098	M486-2	Biopsy cassette	White	£179.00
HIS0099	M486-3	Biopsy cassette	Pink	£179.00
HIS0100	M486-4	Biopsy cassette	Green	£179.00
HIS0101	M486-5	Biopsy cassette	Yellow	£179.00
HIS0102	M486-6	Biopsy cassette	Blue	£179.00
HIS0103	M486-7	Biopsy cassette	Peach	£179.00
HIS0104	M486-8	Biopsy cassette	Tan	£179.00
HIS0105	M486-9	Biopsy cassette	Grey	£179.00
HIS0106	M486-10	Biopsy cassette	Lilac	£179.00
HIS0107	M486-11	Biopsy cassette	Orange	£179.00
HIS0108	M486-12	Biopsy cassette	Aqua	£179.00