

Simport® Scientific Micromesh™ Biopsy Cassettes



This cassette offers 1676 square openings (0.38mm) allowing for a greatly improved fluid exchange and drainage without having to use biopsy pads. Large anterior slots in both cassette and cover ensure that the Micromesh™ biopsy cassette will sink rapidly. A large square compartment measuring 27mm is perfect even for needle biopsies. The cover does not protrude above the cassette, a great space saving feature allowing more cassettes to be stacked in automatic labelling machines and tissue processors.

- Moulded from acetal plastic
- Keep specimens safely submerged in liquid
- Totally resistant to the chemical action of histological solvents
- Anterior writing area is at a 45° angle to make the cassette
- Suitable for use with some automated cassette printers
- Recessed, easy to remove, pre-attached disposable cover
- Available with either one or four compartments



Supplied in packs of 1000 (4 dispenser boxes of 250 pre-assembled cassettes)

Code	Alt Ref	Compartments	Colour	Price
HIS0131	M507-2	1	White	£332.00
HIS0132	M507-3	1	Pink	£332.00
HIS0133	M507-4	1	Green	£332.00
HIS0134	M507-5	1	Yellow	£332.00
HIS0135	M507-6	1	Blue	£332.00
HIS0136	M507-7	1	Peach	£332.00
HIS0137	M507-8	1	Tan	£332.00
HIS0138	M507-9	1	Grey	£332.00
HIS0139	M507-10	1	Lilac	£332.00
HIS0140	M507-11	1	Orange	£332.00
HIS0141	M507-12	1	Aqua	£332.00

Code	Alt Ref	Compartments	Colour	Price
HIS0142	M508-2	4	White	£338.00
HIS0143	M508-3	4	Pink	£338.00
HIS0144	M508-4	4	Green	£338.00
HIS0145	M508-5	4	Yellow	£338.00
HIS0146	M508-6	4	Blue	£338.00
HIS0147	M508-7	4	Peach	£338.00
HIS0148	M508-8	4	Tan	£338.00
HIS0149	M508-9	4	Grey	£338.00
HIS0150	M508-10	4	Lilac	£338.00
HIS0151	M508-11	4	Orange	£338.00
HIS0152	M508-12	4	Aqua	£338.00

Simport® Scientific Microsette® I Biopsy Cassettes



No biopsy pads necessary with these innovative disposable plastic biopsy cassettes. Perfect even with needle biopsies since mesh holes have a diameter of only 0.26mm while still allowing for maximum fluid exchange and drainage.

- Moulded from acetal plastic
- Keep specimens safely submerged in liquid
- Resistant to most histological solvents
- Anterior writing area is at a 45° angle
- Suitable to be used with certain types of cassette labelling instruments
- Over 2000, 0.26mm square openings to maximise fluid exchange
- Lid is pre-attached to base



Supplied in packs of 1000 (4 dispenser boxes of 250 pre-assembled cassettes)

Code	Alt Ref	Compartments	Colour	Price
HIS0109	M502-2	1	White	£360.00
HIS0110	M502-3	1	Pink	£360.00
HIS0111	M502-4	1	Green	£360.00
HIS0112	M502-5	1	Yellow	£360.00
HIS0113	M502-6	1	Blue	£360.00
HIS0114	M502-7	1	Peach	£360.00
HIS0115	M502-8	1	Tan	£360.00
HIS0116	M502-9	1	Grey	£360.00
HIS0117	M502-10	1	Lilac	£360.00
HIS0118	M502-11	1	Orange	£360.00
HIS0119	M502-12	1	Aqua	£360.00

Code	Alt Ref	Compartments	Colour	Price
HIS0120	M503-2	6	White	£360.00
HIS0121	M503-3	6	Pink	£360.00
HIS0122	M503-4	6	Green	£360.00
HIS0123	M503-5	6	Yellow	£360.00
HIS0124	M503-6	6	Blue	£360.00
HIS0125	M503-7	6	Peach	£360.00
HIS0126	M503-8	6	Tan	£360.00
HIS0127	M503-9	6	Grey	£360.00
HIS0128	M503-10	6	Lilac	£360.00
HIS0129	M503-11	6	Orange	£360.00
HIS0130	M503-12	6	Aqua	£360.00