

Prestige Series 2100 Portable Autoclaves



Suitable in nearly all areas of infection control. Designed to the latest international standards for safety and performance and fully CE marked to the latest Medical Devices Directive 93142/EEC (MDD).

- Choice of sterilising cycle
- 126°C for 11 mins
- Two stage overpressure protection
- Operating pressure 1.4 bar
- Supplied complete with a certificate of compliance, performance test certificate, instruction manual and electrical cord set
- AUT4500 Standard Autoclave – complete with small contents plus instrument basket, instrument tray, lifting device
- AUT4502 Extended Autoclave – standard contents plus large instrument basket
- AUT4504 Extended Plus Autoclave – standard contents plus large instrument basket. Supplied with temperature and pressure gauges and thermocouple entry port



Code	Alt Ref	Description	Capacity, L	Max Dims, w x h, mm	Chamber Dims, Ø x d, mm	Max Load, kg	Price
AUT4500	210001	Standard	9	340 x 335	210 x 328	3	£1077.00
AUT4502	210002	Extended	12	340 x 420	210 x 328	4	£1320.00
AUT4504	210004	Extended Plus	12	340 x 420	210 x 328	4	£1450.00

Prestige Classic Media Portable Autoclaves



If you need to sterilise liquids and media or eliminate the risk of cross infection through contaminated waste, then the Classic Media autoclave offers a professional, safe and affordable solution.

- Interlock safety mechanism prevents opening until a safe temperature is reached
- Simple to operate - set the programme with the one touch button, full automatic cycle
- 121°C sterilising temperature
- 1.05 bar operating pressure
- Optional printer for quick and easy data recording
- Special in-built software ensures exact temperatures are reached
- Fully CE marked to the Medical Devices Directive 93/42/EC



Code	Alt Ref	Capacity, L	Max Dims, w x h, mm	Chamber Dims, Ø x d, mm	Max Load, kg	Price
AUT4515	210047UK	9	340 x 335	210 x 228	3	£1240.00
AUT4517	210048UK	12	340 x 420	210 x 318	4	£1490.00

Rodwell Portaclave®™ Portable Autoclaves



Compact and fully portable autoclaves with a robust space saving design. Available with manual control, or with the Vario digital controller which gives variable temperature of 100 to 121°C, variable time of 1 to 99 minutes, and four preset cycles.

- Robust space saving design
- Single phase "plug & play"
- No plumbing needed
- Sodium chloride must not exceed a level greater than 500PPM
- Minimal maintenance
- Grade 316 stainless steel chamber
- Suitable for culture media
- Use regular mains water only, do not use distilled or RO water



Code	Alt Ref	Capacity, L	Dims, w x d x h, mm	Chamber Dims, Ø x d, mm	Controller	Price
AUT2508	ST1528E	15	370 x 390 x 485	280 x 245	Manual	£2350.00
AUT2514	ST1528V	15	370 x 390 x 485	280 x 245	Digital (Vario)	£2950.00
AUT2510	ST2228E	22	370 x 390 x 600	280 x 360	Manual	£2850.00
AUT2516	ST2228V	22	370 x 390 x 600	280 x 360	Digital (Vario)	£3500.00
AUT2512	ST3028E	30	370 x 390 x 730	280 x 490	Manual	£3301.00
AUT2518	ST3028V	30	370 x 390 x 730	280 x 490	Digital (Vario)	£3950.00