

IKA ULTRA-TURRAX® Homogenisers Continued

IKA®

T 65 Dispersing Tools

Code	Alt Ref	Description	Material	Max Temp, °C	Working Volume, mL	Stator/Rotor Ø, mm	Pack	Price
HOM2454	0008005900	S 65 KG - HH - G 65 F	FFPM, WC, AISI 316L	180	2000 - 30000	65/58	Each	£5794.00
HOM2087	0010000643	S 65 KR - G 65 F	FFPM, PTFE, AISI 316L	180	2000 - 30000	65/58	Each	£5262.00
HOM2452	0008005700	S 65 KG - HH - G 65 M	FFPM, WC, AISI 316L	180	2000 - 40000	65/58	Each	£5577.00
HOM2086	0010000642	S 65 KR - G 65 M	FFPM, PTFE, AISI 316L	180	2000 - 40000	65/58	Each	£5048.00
HOM2450	0008005500	S 65 KG - HH - G 65 G	FFPM, WC, AISI 316L	180	2000 - 50000	65/53	Each	£5226.00
HOM2085	0010000641	S 65 KR - G 65 G	FFPM, PTFE, AISI 316L	180	2000 - 50000	65/53	Each	£4695.00



Cole-Parmer® SHM3 Digital Homogeniser

Cole-Parmer®

Can be stand mounted for larger samples, or used as a handheld unit for processing smaller batches. Digital speed control helps increase repeatability, and the soft start motor reduces splashing, keeping both you and your samples safer. Utilise the built in timer to process for seconds to minutes – also features a 30 second countdown for remote applications.

- Powerful and versatile for all your homogenisation needs
- Ideal for tissue homogenisation, cell disruption, emulsions, suspensions, and protein extraction
- Digital control for better repeatability and accuracy
- Process samples from 200µL to 2L
- Vary the 850W motor from 5000 to 25000RPM

Code	Alt Ref	Model	Sample Size, mL	Speed Range, RPM	Dims, h x Ø, mm	Weight, kg	Price
HOM2508	04720-00	SHM3	0.2 - 2000	5000 - 25000	70 x 216	1.6	£2199.00

Rotor Probe

Code	Alt Ref	Capacity, mL	Dims, l x Ø, mm	Price
HOM2510	04720-03	0.03 - 5	75 x 5	£919.00

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty