

Corning® CoolRacks™

CORNING

Thermo-conductive modules which reduce variability which originates from tubes or wells placed directly into ice, dry ice, alcohol baths, water baths, and other temperature sources. Place the thermo-conductive modules directly onto a temperature source between -196°C to $>100^{\circ}\text{C}$ and it will rapidly adapt to that temperature. Both CoolRack™ and CoolSink® modules ensure $\pm 0.1^{\circ}\text{C}$ temperature uniformity across all tubes and plates when cooling, snap-freezing, heating, or thawing. All thermo-conductive modules may be autoclaved, high heat sterilised, or decontaminated with bleach, alcohol, or other disinfectants or lab detergents.

Applications include:

- Cooling reagents such as ECMs, like Corning® Matrigel® matrix
- Restriction enzymes
- dNTPs and antibodies
- Alcohol-free dry ice snap-freezing of tissue



A



A



A

Code	Description	Capacity	Colour	Price
432034	CoolRack M6	6 x 1.5/2mL microcentrifuge tubes	Grey	£69.00
432035	CoolRack M6	6 x 1.5/2mL microcentrifuge tubes	Green	£69.00
432036	CoolRack M6	6 x 1.5/2mL microcentrifuge tubes	Orange	£69.00
432037	CoolRack M15	15 x 1.5/2mL microcentrifuge tubes	Grey	£162.00
432038	CoolRack M15	15 x 1.5/2mL microcentrifuge tubes	Green	£162.00
432039	CoolRack M15	15 x 1.5/2mL microcentrifuge tubes	Orange	£162.00
432040	CoolRack XT M24	24 x 1.5/2mL microcentrifuge tubes	Grey	£280.00
432041	CoolRack M30	30 x 1.5/2mL microcentrifuge tubes	Grey	£282.00
432042	CoolRack M30	30 x 1.5/2mL microcentrifuge tubes	Green	£282.00
432043	CoolRack M30	30 x 1.5/2mL microcentrifuge tubes	Orange	£282.00
432044	CoolRack M90	90 x 1.5/2mL microcentrifuge tubes	Grey	£426.00
432045	CoolRack M96 with A-H and 1-12 row and column indexing	96 x 1.5/2mL microcentrifuge tubes	Grey	£562.00
432046	CoolRack M30-PF	30 x 500µL conical microcentrifuge tubes	Grey	£333.00
432047	CoolRack M15-PF	15 x 1.5mL conical microcentrifuge tubes	Grey	£222.00
432048	CoolRack M30-PF	30 x 1.5mL conical microcentrifuge tubes	Grey	£333.00
432049	CoolRack CF15	15 x cryovial/FACS tubes	Grey	£162.00
432050	CoolRack XT CFT24 with "gripping" wells for one hand vial opening/closing	24 x cryovial/FACS tubes	Grey	£319.00
432052	CoolRack CFT30 with "gripping" wells for one hand vial opening/closing	30 x cryovial/FACS tubes	Grey	£333.00
432051	CoolRack CF45	45 x cryovial/FACS tubes	Grey	£301.00
432053	CoolRack XT PCR96	12 x strip wells	Grey	£335.00
432054	CoolRack XT-M-PCR	12 x 1.5/2mL microcentrifuge tubes and 6 x strip wells	Grey	£327.00
432055	CoolRack XT PCR384	384 well plate	Grey	£372.00
432056	CoolRack 96x0.5mL	96 x 0.5mL 2D tubes	Grey	£413.00
432057	CoolRack 96x1mL	96 x 1.4mL 2D tubes	Grey	£378.00
432058	CoolRack SV2	12 x 5mL injectable ampules	Grey	£347.00
432059	CoolRack SV10	12 x 10mL injectable ampules	Grey	£347.00
432060	CoolRack XT 5mL	12 x 5mL centrifuge tubes	Grey	£282.00
432061	CoolRack 15mL	9 x 15mL conical centrifuge tube	Grey	£296.00
432062	CoolRack 50mL	4 x 50mL conical centrifuge tube	Grey	£310.00
432063	CoolRack 250mL	1 x 250mL conical centrifuge tube	Grey	£335.00