

Corning® CoolRacks™ Continued



Code	Description	Capacity	Colour	Price
432064	CoolRack 250mL	1 x 250mL Easy Grip storage bottle	Grey	£335.00
432065	CoolRack VS13	9 x 13 x 75mm blood tubes	Grey	£309.00
432066	CoolRack V16	9 x 16 x 100mm blood tubes	Grey	£309.00
432067	CoolRack V13	9 x 13 x 100mm blood tubes	Grey	£309.00
432069	CoolRack LV lightweight insulated module	12 x 13/16mm diameter tubes	Grey	£292.00
432068	CoolRack L lightweight insulated module	12 x 15mL centrifuge tubes	Grey	£292.00

Corning® CoolBox™ Systems



Corning® CoolBox™ Systems are bench top cooling systems that provide sample cooling or freezing without ice, electricity, or batteries. Versatile and accommodating to a variety of sample formats and temperatures, the CoolBox™ modular design enables the use of both Corning® CoolRack® and Corning® CoolSink® thermo-conductive sample modules to hold microcentrifuge tubes, cryogenic vials, PCR tubes, or assay plates at a stable temperature while you work.

Corning® XT Starter is a low profile, open-platform cooler that accommodates most Corning® CoolRack® and Corning® CoolSink® modules, making it a versatile and flexible tool for a variety of applications. XT Starter keeps samples cold (0.5 to 4.0°C) for up to four hours. The 1 to 8°C temperature indicator provides visual assurance of proper temperature.

Corning® CoolBox™ XT systems are ideal as a benchtop cooling workstation. This system includes a base, collar, lid, and one or two Corning® XT Cooling Cores which can be used for sample cooling and accommodates one tube or plate sample module.



Code	Description	Colour	Supplied with	Price
432014	XT starter, complete	Purple	XT holder, XT cooling core	£221.00
432015	XT starter, holder only	Purple	XT holder	£68.00
432021	CoolBox XT	Purple	Holder, cooling core	£537.00
432022	CoolBox XT	Green	Holder, cooling core	£537.00
432023	CoolBox XT	Orange	Holder, cooling core	£537.00
432024	CoolBox XT	Pink	Holder, cooling core	£537.00
432025	CoolBox XT	Purple	Holder, 2 x cooling cores	£825.00
432026	CoolBox XT	Green	Holder, 2 x cooling cores	£825.00
432027	CoolBox XT	Orange	Holder, 2 x cooling cores	£825.00
432028	CoolBox XT	Pink	Holder, 2 x cooling cores	£825.00

Accessories

Code	Description	Colour	Price
432081	XT cooling core for CoolBox XT Workstations for 0.5 to 4°C	—	£178.00
432082	XT freezing core for CoolBox XT Workstations for below 0°C	—	£244.00
432083	Extension collar for CoolBox XT	Purple	£94.00
432084	Extension collar for CoolBox XT	Green	£94.00
432085	Extension collar for CoolBox XT	Orange	£94.00
432086	Extension collar for CoolBox XT	Pink	£94.00
432088	Extension collar for CoolBox 2XT	Green	£117.00
432089	Extension collar for CoolBox 2XT	Orange	£117.00
432090	Extension collar for CoolBox 2XT	Pink	£117.00
432091	1 to 8°C adhesive strip LCD temperature display	—	£50.00
432092	Elastic sleeves for CoolRack modules	—	£28.00