

Genlab Multi-Purpose Benchtop Incubators Continued



Aluminium Coated Mild Steel Chamber

Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
INC2008	MINI/6/TDIG	6	1	1	350 x 330 x 460	230 x 190 x 150	£891.00
INC2010	MINI/18/TDIG	18	1	1	380 x 410 x 550	260 x 260 x 270	£946.00
INC2012	MINI/30/TDIG	30	2	2	480 x 490 x 510	360 x 350 x 240	£1002.00
INC2014	MINI/40/TDIG	40	2	3	480 x 490 x 540	360 x 350 x 320	£1085.00
INC2016	MINI/50/TDIG	53	2	3	630 x 470 x 660	490 x 330 x 330	£1157.00
INC2018	MINI/75/TDIG	73	2	3	630 x 600 x 660	490 x 450 x 450	£1290.00
INC2020	MINI/100/TDIG	99	3	4	630 x 600 x 780	490 x 450 x 450	£1423.00
INC2022	MINI/125/TDIG	121	3	5	630 x 600 x 880	490 x 450 x 550	£1557.00
INC2024	MINI/150/TDIG	143	3	5	630 x 680 x 880	490 x 530 x 550	£1713.00
INC2026	MINI/200/F/TDIG	180	4	7	630 x 840 x 1080	490 x 490 x 750	£2185.00

Stainless Steel Chamber

Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
INC2028	MINI/6/TDIG/SS	6	1	1	350 x 330 x 460	230 x 190 x 150	£940.00
INC2030	MINI/18/TDIG/SS	18	1	1	380 x 410 x 550	260 x 270 x 260	£1002.00
INC2032	MINI/30/TDIG/SS	30	2	2	480 x 490 x 510	360 x 350 x 240	£1067.00
INC2034	MINI/40/TDIG/SS	40	2	3	480 x 490 x 540	360 x 350 x 320	£1168.00
INC2036	MINI/50/TDIG/SS	53	2	3	63 x 47 x 66	490 x 330 x 330	£1250.00
INC2038	MINI/75/TDIG/SS	73	2	3	630 x 600 x 660	490 x 450 x 330	£1396.00
INC2040	MINI/100/TDIG/SS	99	3	4	630 x 600 x 780	490 x 450 x 450	£1541.00
INC2042	MINI/125/TDIG/SS	121	3	5	630 x 600 x 880	490 x 450 x 550	£1696.00
INC2044	MINI/150/TDIG/SS	143	3	5	630 x 680 x 880	490 x 530 x 550	£1670.00
INC2046	MINI/200/F/TDIG/SS	180	4	7	630 x 840 x 1080	490 x 530 x 550	£2418.00



Genlab Prime Incubators



The flagship brand of Genlab Incubators. They are built to the highest standards and packed with many features and enhancements. As standard, the Prime range offers a mirror polished stainless steel interior, increased accuracies and an extended temperature range of 100°C. It also features the latest touch screen control system which offers intuitive control designed for laboratory incubators.

- Mirror polished stainless steel internal chamber
- Temperature range: ambient +5°C to 100°C
- Fluctuation: $\pm 0.15^\circ\text{C}$ @ 37°C
- On screen historical trending (48 hours)
- 8 stage profile control
- Countdown timer
- 24/7 real time clock operation
- Increased ramp up/recovery times
- Manual or automatic overheat reset
- High accuracy Pt100A duplex sensors $< 0.6^\circ\text{C}$
- Fan assisted circulation
- Excellent stability $< \pm 0.5^\circ\text{C}$
- Low chamber uniformity

Code	Alt Ref	Capacity, L	No. of Shelves	Shelf Positions	Ext Dims, w x d x h, mm	Price
INC2048	PRI/30/TDIG	30	2	2	500 x 620 x 570	£1375.00
INC2058	PRI/50/TDIG	50	2	3	630 x 590 x 660	£1552.00
INC2224	PRI/75/TDIG	75	2	3	630 x 720 x 660	£1707.00
INC2226	PRI/100/TDIG	100	3	4	630 x 720 x 780	£1807.00
INC2228	PRI/125/TDIG	125	3	5	630 x 720 x 880	£1963.00
INC2230	PRI/150/TDIG	150	3	5	630 x 800 x 880	£2124.00
INC2232	PRI/200/TDIG	200	4	7	630 x 800 x 1080	£2529.00