

OHAUS Adventurer® Precision Balances Continued



Code	Alt Ref	Model	Capacity, g	Readability, g	Calibration	Pan Dims, mm	Dims, w x d x h, mm	Price
BAL0016	30127961	AX223/E	220	0.001	External	Ø 130	354 x 230 x 340	£1335.00
BAL0015	30127960	AX223	220	0.001	Internal	Ø 130	354 x 230 x 340	£1615.00
BAL0018	30127963	AX423/E	420	0.001	External	Ø 130	354 x 230 x 340	£1495.00
BAL0017	30127962	AX423	420	0.001	Internal	Ø 130	354 x 230 x 340	£1780.00
BAL0020	30122615	AX523/E	520	0.001	External	Ø 130	354 x 230 x 340	£1840.00
BAL0019	30122614	AX523	520	0.001	Internal	Ø 130	354 x 230 x 340	£2120.00
BAL0056	30206080	AX422/E	420	0.01	External	175 x 195	354 x 230 x 100	£1015.00
BAL0054	30206079	AX422	420	0.01	Internal	175 x 195	354 x 230 x 100	£1305.00
BAL0060	30206082	AX822/E	820	0.01	External	175 x 195	354 x 230 x 100	£1165.00
BAL0058	30206081	AX822	820	0.01	Internal	175 x 195	354 x 230 x 100	£1445.00
BAL0024	30127967	AX1502/E	1520	0.01	External	175 x 195	354 x 230 x 100	£1315.00
BAL0023	30127966	AX1502	1520	0.01	Internal	175 x 195	354 x 230 x 100	£1595.00
BAL0026	30127969	AX2202/E	2200	0.01	External	175 x 195	354 x 230 x 100	£1345.00
BAL0025	30127968	AX2202	2200	0.01	Internal	175 x 195	354 x 230 x 100	£1625.00
BAL0028	30122617	AX4202/E	4200	0.01	External	175 x 195	354 x 230 x 100	£1475.00
BAL0027	30122616	AX4202	4200	0.01	Internal	175 x 195	354 x 230 x 100	£1760.00
BAL0029	30122618	AX5202	5200	0.01	Internal	175 x 195	354 x 230 x 100	£2120.00
BAL0064	30206084	AX2201/E	2200	0.1	External	175 x 195	354 x 230 x 100	£970.00
BAL0062	30206083	AX2201	2200	0.1	Internal	175 x 195	354 x 230 x 100	£1245.00
BAL0068	30122620	AX8201/E	8200	0.1	External	175 x 195	354 x 230 x 100	£1290.00
BAL0066	30122619	AX8201	8200	0.1	Internal	175 x 195	354 x 230 x 100	£1575.00

Accessories

Code	Alt Ref	Description	Price
BAL5162	80251395	Auxiliary display, TAD7	£48.10
BAL6434	12106964	RS232 cable PC 9 pin, female to female	£75.00
BAL3954	80850000	Anti-theft device	£54.60
BAL3959	83020221	In use display cover	£18.30

A&D FXi Series Precision Balances



The FXi series offers high performance and quality combined with accessibility, bringing to you a cost effective solution for a diverse array of applications.

- 2nd generation super hybrid sensor
- Die cast and stainless steel body
- High speed response time for stability in 1 second
- Blue vacuum fluorescent display
- Large round weighing pan
- Arithmetical statistics function
- Standard RS-232C interface
- USB and LAN options
- Compact size: 263 x 193 x 85mm



Code	Alt Ref	Capacity, g	Readability, g	Pan Dims, mm	Price
BAL4600	FX-120i	122	0.001	Ø 130	£566.00
BAL4601	FX200i	220	0.001	Ø 130	£609.00
BAL4602	FX-300i	320	0.001	Ø 130	£646.00
BAL4603	FX-500i	520	0.001	Ø 130	£730.00
BAL4604	FX-1200i	1220	0.01	Ø 150	£566.00
BAL4605	FX-2000i	2200	0.01	Ø 150	£609.00
BAL4606	FX-3000i	3200	0.01	Ø 150	£646.00
BAL4607	FX-5000i	5200	0.01	Ø 150	£725.00