

LMS Series 1A Cooled Incubators

LMS

Reliable and accurate temperature-controlled cabinets, operating between -10 to +50°C. With 4 model sizes to choose from these versatile cabinets will suit most workloads ranging from 76 to 272L capacity.

- Digital set/display
- Fan-assisted circulation
- Hermetically-sealed refrigeration system
- Self-evaporating condensation drip tray
- Preset overheat protection
- For 220/240V, single phase supply
- White pre-coated aluminium interior
- A full range of upgrades including fruit fly modification are also available to suit your requirements, please enquire for details



Single Set Point Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5300	80	76	2	520 x 680 x 780	£2594.00
INC5302	120	120	3	520 x 680 x 1080	£2965.00
INC5304	201	201	4	520 x 680 x 1430	£3200.00
INC5306	280	272	5	520 x 680 x 1750	£3336.00

Automatic Dual Temperature Control

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5308	80A	76	2	520 x 680 x 780	£3004.00
INC5310	120A	120	3	520 x 680 x 1080	£3375.00
INC5312	201A	201	4	520 x 680 x 1430	£3610.00
INC5314	280A	272	5	520 x 680 x 1750	£3746.00

Fully Programmable Controller (Multiple Temperatures/Time Frames)

Code	Alt Ref	Capacity, L	No. of Shelves	Ext Dims, w x d x h, mm	Price
INC5316	80NP	76	2	520 x 680 x 780	£2944.00
INC5318	120NP	120	3	520 x 680 x 1080	£3315.00
INC5320	201NP	201	4	520 x 680 x 1430	£3550.00
INC5322	280NP	272	5	520 x 680 x 1750	£3686.00

LMS Series 2 Cooled Incubators

LMS

Offering reliable and accurate temperature-control, operating between -10 to +50°C. The range comprises of 4 cabinet sizes to suit most workloads, ranging from 135 to 420L capacity. With extra width these cabinets offer better access to the internal chamber.

- Digital set/display
- Fan-assisted circulation
- Hermetically-sealed refrigeration system
- Self-evaporating condensation drip tray
- Preset overheat protection
- For 220/240V, single phase supply
- ABS plastic interior
- A full range of upgrades including fruit fly modification are also available to suit your requirements, please enquire for details

