

Memmert ICP Compressor Cooled Incubators Continued

Code	Model	Volume, L	Int Dims*, w x d x h, mm	Ext Dims**, w x d x h, mm	Price
INC8120	ICP110	108	560 x 400 x 480	745 x 584 x 1233	£5370.00
INC8122	ICP260	256	640 x 500 x 800	824 x 684 x 1552	£8008.00
INC8124	ICP450	449	1040 x 600 x 720	1224 x 784 x 1467	£9456.00
INC8126	ICP750	749	1040 x 600 x 1200	1224 x 784 x 1950	£11315.00

* Internal dimensions: depth less 33mm for fan

** External dimensions: depth without door handle, depth of handle 56mm, height with castors

Accessories

Further options and accessories are available on request. Additional equipment must be specified when ordering.

Code	Alt Ref	Description	For use with	Price
CHA4022	E20165	Stainless steel grid, electropolished	ICP110	£134.00
CHA4052	E28891	Stainless steel grid, electropolished	ICP260	£160.00
CHA4054	E20182	Stainless steel grid, electropolished	ICP450, ICP750	£246.00
CHA4062	B29725	Perforated stainless steel shelf	ICP260	£171.00
CHA4064	B00328	Perforated stainless steel shelf	ICP450, ICP750	£261.00



Memmert ICPeco Compressor Cooled Incubators

This environmentally-friendly incubator is cooled with climate-neutral CO₂. Thanks to the excellent thermodynamic properties of the refrigerant CO₂ (R744) and the finely adjusted control technology, the ICPeco cooled incubator is both powerful and high-precision. Without critical temperature overshoots, it keeps the temperatures exactly at the setpoint.

- Temperature range: -12 to +60°C
- 3 model sizes
- TwinDISPLAY - adaptive multifunctional digital PID-microprocessor controller with 2 high-definition TFT-colour displays
- Double doors standard for all models: Prevention of contamination and drops in temperature, and at the same time an optimum view of the sensitive load through wide-area interior glass doors
- Very good air circulation in the working chamber through sophisticated ventilation technology
- Intelligent defrosting function

For further information on TwinDISPLAY please see page 529

Code	Model	Capacity, L	Int Dims*, w x d x h, mm	Ext Dims**, w x d x h, mm	Price
INC2400	ICP260eco	256	640 x 500 x 800	824 x 684 x 1552	£9253.00
INC2402	ICP450eco	449	1040 x 600 x 720	1224 x 784 x 1467	£11348.00
INC2404	ICP750eco	749	1040 x 600 x 1200	1224 x 784 x 1950	£13173.00

* Internal dimensions: depth less 33mm for fan

** External dimensions: depth without door handle, depth of handle 56mm