



A&D FXi-WP Series Water/Dust Proof Precision Balances



The FXi-WP series has all the benefits of the FXi series, with the additional feature of being IP65 dust and waterproof making them ideal for hostile environments. This enables the measurement of powder or liquid at 1mg resolution, and eliminates any worries about damage to the balance.

- IP-65 dust & water protection, ideal for food labs, teaching labs etc.
- 2nd generation super hybrid sensor
- Die cast and stainless steel body
- High speed response time for stability in 1 second
- Blue vacuum fluorescent display
- Large round weighing pan
- Arithmetical statistics function
- Optional IP-65 RS-232C cable
- USB and LAN options (IP-65 not available for these options)
- Compact size: 263 x 193 x 85mm

| Code | Alt Ref | Capacity, g | Readability, g | Pan Dims, mm | Price |
|---------|-------------|-------------|----------------|--------------|----------|
| BAL2090 | FX-120I-WP | 122 | 0.001 | Ø 130 | £931.00 |
| BAL2092 | FX-300I-WP | 320 | 0.001 | Ø 130 | £1028.00 |
| BAL2094 | FX-1200I-WP | 1220 | 0.01 | Ø 150 | £931.00 |
| BAL2096 | FX-3000I-WP | 3200 | 0.01 | Ø 150 | £1028.00 |

Options and Accessories

| Code | Alt Ref | Description | Price |
|---------|-----------|---|----------|
| BAL4650 | Fxi-OP-02 | Quick USB interface with cable | £98.00 |
| BAL4651 | Fxi-OP-08 | LAN ethernet interface | £490.00 |
| BAL4652 | Fxi-OP-09 | Internal rechargeable NiMH battery pack | £358.00 |
| BAL4653 | Fxi-OP-10 | Small breeze break with doors | £65.00 |
| BAL4656 | AD-1672 | Tabletop breeze break | £1942.00 |



A&D GX-A/GF-A Series Precision Balances



A range of precision top pan balances from A&D Instruments with innovative features especially suited for pharmaceutical and life science applications.

- Patented electronically controlled load (ECL) mode for automatic determination of minimum weight and standard deviation for GLP and FDA applications
- Flow Rate Display (FRD) mode for displaying flow rate (e.g. in mL/min) from a pump or dosing system, high and low level warning can be set using comparator mode
- Minimum weight alert on display when weighing below calculated or input minimum weight
- 10 user level passwords, 1 supervisor password
- Quick and bidirectional USB options with built in USB port
- Impact Shock Detection (ISD) warns user when balance is being potentially damaged due to impact load shocks. Up to 50 impacts stored in balances memory with date and time
- Fast stabilisation, typically 1 second
- Automatic self calibration (GX models only)
- Under hook for below balance weighing e.g. density determination
- Free downloadable data capture and control software

| Code | Alt Ref | Capacity, g | Readability, g | Pan Dims, mm | Calibration | Price |
|---------|---------|-------------|----------------|--------------|-------------|----------|
| BAL1378 | GF123A | 120 | 0.001 | 128 x 128 | External | £1028.00 |
| BAL1382 | GF203A | 220 | 0.001 | 128 x 128 | External | £1136.00 |
| BAL1380 | GX203A | 220 | 0.001 | 128 x 128 | Internal | £1456.00 |
| BAL1386 | GF303A | 320 | 0.001 | 128 x 128 | External | £1217.00 |
| BAL1384 | GX303A | 320 | 0.001 | 128 x 128 | Internal | £1536.00 |
| BAL1390 | GF403A | 420 | 0.001 | 128 x 128 | External | £1271.00 |