


SLS Select Loops

- For inoculation of nutrient media
- Available aseptically manufactured or sterile by irradiation



Supplied in packs of 1000 (50 packs of 20)


Code	Volume, µL	Description	Colour	Sterility	Price
SLS2008	1	Flexible	Blue	Non sterile	£44.00
SLS2014	1	Firm	Clear	Sterile	£82.00
SLS2010	10	Flexible	Blue	Non sterile	£42.90
SLS2012	10	Firm	Green	Sterile	£82.00

SLS Select Spreaders



- Polystyrene spreaders
- Available aseptically manufactured or sterile by irradiation

Code	Description	Colour	Sterility	Inner Pack	Pack	Price
SLS2098	L-shaped	Blue	Non sterile	5	500	£46.20
SLS2102	L-shaped	Blue	Non sterile	20	500	£57.10
SLS2104	L-shaped	Clear	Sterile	1	1000	£202.00
SLS1999	Wedge-shaped	Blue	Non sterile	1	250	£43.80




Greiner Spreaders

- Polystyrene L-shaped spreaders
- Supplied sterile (SAL 10⁻³)
- Dims, l x w: 152 x 38mm

Supplied in packs of 500


Code	Inner Pack	Price
G730190	20	£21.20
G730191	5	£35.40



Thermo Scientific™ Sterilin™ Loops and Spreaders



- For inoculation, dilution streaking and spreading
- Eliminate the need for flaming and the subsequent risk of aerosols
- Choice of rigid or flexible loops to suit different applications and user preferences
- Hexagonal shaft improves grip and assists orientation
- Free of lubricants, oils and electrostatic charges to facilitate consistent wetting and complete liquid transfer
- No rough edges means smooth, problem-free plating and streaking of cultures
- Sterilised by gamma irradiation



Code	Alt Ref	Volume, µL	Description	Colour	Inner Pack	Pack	Price
IN02000	SL1S	1	Loops, soft	Light green	10	800	£105.00
IN02001	SL1H	10	Loops, hard	Dark green	10	800	£105.00
IN02004	QL1	10	Loops, with sphere	Green	10	1000	£154.00
IN02002	SL10S	5	Loops, soft	Light blue	20	800	£105.00
IN02003	SL10H	10	Loops, hard	Dark blue	10	800	£111.00
IN02006	QL10	10	Loops, with sphere	Blue	20	1000	£162.00
X7360	SPCS01	—	Spreader, L-shaped	Blue	1	500	£266.00