



OHAUS Navigator™ Portable Balances



Rugged, multi-purpose balances for a wide range of general weighing applications, designed for accuracy and repeatability in essential weighing applications. Improved productivity with fast stabilisation time and mechanical overload protection safeguards the balance increasing durability.

- 4 applications
- Backlit liquid crystal display (LCD)
- AC adaptor (included) or 4 x C batteries (not included)
- Communication RS232, USB, or ethernet (available as accessories)
- Constructed of ABS housing, stainless steel pan, transportation lock, slip resistant and adjustable feet, levelling bubble
- Design features include menu and calibration lock switch, mechanical and software overload/underload protection, stability indicator, low battery indicator, auto shut off, user selectable print options, user selectable communication settings

Code	Alt Ref	Model	Capacity, g	Readability, g	Pan Dims, mm	Price
BAL1280	30615581	NV123	120	0.001	Ø 93	£236.00
BAL1282	30615582	NV223	220	0.001	Ø 93	£264.00
BAL1284	30615583	NV323	320	0.001	Ø 93	£278.00
BAL1426	30456382	NV222	220	0.01	Ø 145	£239.00
BAL1428	30456383	NV422	420	0.01	Ø 145	£246.00
BAL1430	30456384	NV622	620	0.01	Ø 145	£251.00
BAL1286	30615584	NV1202	1200	0.01	190 x 144	£256.00
BAL1288	30615585	NV2202	2200	0.01	190 x 144	£268.00
BAL1290	30615586	NV3202	3200	0.01	190 x 144	£285.00
BAL1432	30456385	NV221	220	0.1	190 x 144	£144.00
BAL1434	30456386	NV621	620	0.1	190 x 144	£170.50
BAL1436	30456387	NV1201	1200	0.1	190 x 144	£196.50
BAL1438	30456388	NC2201	2200	0.1	190 x 144	£207.00
BAL1440	30456389	NVT2201	2200	0.1	230 x 174	£219.00
BAL1442	30456390	NVT4201	4200	0.1	230 x 174	£229.00
BAL1444	30456391	NVT6201	6200	0.1	230 x 174	£241.00
BAL1292	30615587	NVT10201	10200	0.1	230 x 174	£251.00
BAL1450	30456392	NVT2200	2200	1	230 x 174	£152.50
BAL1452	30456393	NVT6200	6200	1	230 x 174	£175.50
BAL1446	30456394	NVT12000	12000	1	230 x 174	£233.00
BAL1448	30456395	NVT22000	22000	1	230 x 174	£256.00

Accessories

Code	Alt Ref	Description	Price
BAL5528	83032107	RS232 interface kit	£38.70
BAL5530	83032108	USB interface kit	£38.70

Tough samples?

Check out our range of mills and grinders.

See pages 502–511

