

Laboratory Jacks

- Top and bottom plates manufactured from aluminium, giving exceptional rigid and robust support without excessive weight
- Blue platforms
- Stable and free from rocking movements

Code	Plate Dims, w x d, mm	Max Height, mm	Max Load, kg	Price
JAC0998	50 x 40	78	20	£140.00
JAC1000	100 x 100	140	30	£132.00
JAC1004	150 x 150	290	60	£175.00
JAC1008	200 x 200	290	80	£220.00
JAC1012	200 x 230	290	80	£267.00
JAC1013	320 x 260	470	90	£597.00



Extension Plate

- Anodised aluminium
- Complete with stainless steel fixing screws

Code	Description	For use with	Price
JAC1020	300mm extension plate	JAC1008 and JAC1012 only	£182.00

Laboratory Jack

- Screw operated scissor action mechanism
- Strong, stable support for laboratory equipment
- Finished in blue colour

Code	Max Load, kg	Height, mm	Plate Dims w x l, mm	Price
CH0642	25	60 - 250	150 x 135	£58.00



Take a closer look

at our microscope range!

Including monocular, binocular and trinocular models.

See pages 623–649



SCIENTIFIC
LABORATORY
SUPPLIES

