

Honey Pot Jars

- Polyethylene terephthalate (PET) material
- Wide mouth storage containers
- With secure, screw top lid
- Ideal for storing dry powders and other solid materials



Code	Description	Capacity, mL	Label	Sterility	Pack	Price
CON7773	Round	300	No	Sterile	60	£178.00
CON7777	Square	375	Yes	Non sterile	100	£72.00
CON7775	Square	375	No	Sterile	100	£75.00

Kilner Type Jars

- Wide mouth
- Clear glass with lacquered metal plate with flowed-in plastic seal
- Separate metal screw collar

Code	Capacity, mL	Price
JAR1000	500	£6.00
JAR1002	1000	£6.00



Azlon® PP Jugs with Printed Graduations



- Tapered jugs in high clarity polypropylene
- Ergonomically designed handle with thumb grip for easy carrying
- Suitable for food contact applications in accordance with EU and U.S FDA regulations
- Crisp blue printed graduations
- Autoclavable at 121°C



Code	Alt Ref	Capacity, mL	Graduations, mL	Dims, Ø x h, mm	Pack	Price
JUG1096	JPB0050P	50	10	45 x 68	5	£7.70
JUG1098	JPB0100P	100	10	61 x 72	5	£8.60
JUG1100	JPB0250P	250	50	70 x 120	Each	£2.80
JUG1102	JPB0500P	500	100	90 x 140	Each	£3.10
JUG1104	JPB1000P	1000	100	104 x 180	Each	£4.00
JUG1106	JPB2000P	2000	500	150 x 214	Each	£7.80
JUG1108	JPB3000P	3000	500	152 x 240	Each	£9.70
JUG1110	JPB5000P	5000	500	200 x 500	Each	£15.90