

## SLS Select Cryogenic Labels



- Self adhesive labels designed specifically for low temperature use down to  $-196^{\circ}\text{C}$  (apply at  $5^{\circ}\text{C}$  or more)
- Contain a unique adhesive formulation that guarantees extra staying power when applying to polycarbonate, polypropylene, glass and most metals
- White



### Dots

| Code    | Ø, mm | Pack         | Price  |
|---------|-------|--------------|--------|
| LAB2210 | 9.5   | Reel of 1000 | £28.60 |
| LAB2212 | 9.5   | Reel of 5000 | £67.00 |

### Strips

| Code    | Dims, w x l, mm | Pack                           | Price  |
|---------|-----------------|--------------------------------|--------|
| LAB2200 | 33 x 13         | Reel of 1000                   | £30.50 |
| LAB2230 | 33 x 13         | 1700 (20 sheets of 85 labels)  | £60.00 |
| LAB2232 | 24 x 13         | 2380 (20 sheets of 119 labels) | £60.00 |
| LAB2234 | 67 x 25         | 600 (20 sheets of 30 labels)   | £58.00 |

## Hi-Low Freezer Labels

- Self adhesive vinyl labels
- Can be used with most laser and inkjet printers for clear and quick bulk label creation
- Smear resistant
- Offer a high degree of chemical resistance, humidity and UV resistance
- Can be used on polypropylene and polycarbonate surfaces
- Will not shrink if using at the higher temperature band
- Suitable for use between  $-85$  and  $100^{\circ}\text{C}$ , can also be applied at  $-85^{\circ}\text{C}$
- White



| Code    | Size, mm | Pack                           | Price  |
|---------|----------|--------------------------------|--------|
| LAB2240 | 33 x 13  | 1980 (22 sheets of 90 labels)  | £36.00 |
| LAB2242 | 25 x 13  | 2016 (16 sheets of 126 labels) | £32.40 |

## Marker Pens

- Fine point
- Writes on cold or wet surfaces
- Quick drying

### Supplied in packs of 10

| Code    | Colour | Price  |
|---------|--------|--------|
| LAB2900 | Black  | £54.00 |
| LAB2902 | Blue   | £63.00 |
| LAB2904 | Green  | £63.00 |
| LAB2906 | Red    | £57.00 |

