



Edding Laboratory Markers

- Quick-drying, wear-resistant and waterproof and lightfast
- Permanent on plastics, removable from glass with acetone, ethanol, isopropanol and methyl ethyl ketone
- Extreme resistance to abrasion, high temperatures on aluminium up to 500°C and cold temperatures down to -183°C in liquid nitrogen
- Low odour ink with no added toluene/xylene
- Black

Supplied in packs of 10

Code	Nib, mm	Price
E8014	1	£33.50
E8015F	0.75	£33.50



Brady BMP®21-LAB Label Printers



An all-in-one labeller for research, academic and clinical labs. With label materials designed specifically for harsh environments, the BMP21-LAB printer can quickly and easily create legible labels for flat, curved or highly textured surfaces that stay stuck for years even when exposed to extreme temperatures and chemicals

- Label materials available for cryo vials and slides
- Runs on 6 x AA, batteries; rechargeable Li-Ion battery available as an option (inc. in the starter kit)
- Can print commonly used laboratory symbols with intuitive menu structure
- Vial label function automatically cuts the tape according to the size of the vial
- Automated serialisation, time and date stamp
- 3 line LCD display
- Smart Cell Technology for automatic label setup
- Automatic label formatting for commonly used label types in laboratory
- For continuous labels
- Monocolour
- Resolution: 203dpi

Code	Alt Ref	Description	Supplied with	Price
LAB6016	139538	Blister pack	Quick start guide and M21-750-7425 label cartridge	£120.00
LAB6029	134475	Starter kit	Li-ion battery, power cord, quick start guide and M21-750-7425 label cartridge	£184.00

Accessories

Code	Alt Ref	Description	Price
LAB6017	110417	Power adapter	£55.00
LAB6019	110891	Magnet accessory	£16.06



Brady BBP®12 Label Printers



The BBP®12 packs a lot of punch into one compact tool with its high-speed printing, an extended list of labelling options and a simple, user-friendly design. It gives you the versatility to print all your basic label types from one machine - saving costs, reducing downtime and simplifying your labelling.

- Wide variety of labels available, including external-feed bulk rolls and custom parts, for general identification, wire marking, asset tracking and more
- Specifically designed small-core media rolls and ribbons provide fast installation and reduce change-over time
- Centering ribbon spindle virtually eliminates ribbon wrinkle