

Folding Magnifiers

Folding pocket magnifier with a plano-convex lens in an plastic lens holder with aluminium anodised cover. Lens Ø: 17mm.

Code	Magnification	Price
MAG2030	x8	£7.44
MAG2032	x10	£9.45



INTEGRA Biosciences MEDIACLAVE Media Preparators

INTEGRA

Allows the rapid and gentle sterilisation of 1 – 30L of culture medium. Precise controlling and monitoring of process parameters during the sterilisation process guarantee media of a consistently high quality. The intuitive interface and the simple programming make it very easy, even for untrained users, to operate. Can be conveniently connected to the MEDIAJET, allowing the automated and continuous dispensing of media into Petri dishes at the touch of a button.

- Flow heater system for fast and gentle media preparation
- Powerful bidirectional stirrer to guarantee homogeneous media
- Controlled deaeration to guarantee complete sterilisation
- *In situ* sterilisation of the dispense port to ensure sterility of the dispensed media
- Indication of service intervals to avoid unscheduled downtime
- Intuitive multilingual user interface for easy and convenient use
- TÜV certified for compliance to international safety standards
- Integrated web server for online documentation of the process
- Electronically signed log files according to FDA and EU guidelines
- Supplied with cuvette, lid seal, magnetic stirrer, decanting tubing, and fitting for dispensing tubing
- Models with printer also supplied with integrated printer and printer splash guard



Code	Model	Capacity, L	Dims, w x d x h, mm	Weight, kg	Power Supply	Price
I136025	MC 10	10	550 x 640 x 480	57	200-240V, 50/60Hz, 16A, no plug	£22888.00
I136020	MC 10 with printer	10	550 x 640 x 480	57	200-240V, 50/60Hz, 16A, no plug	£24495.00
I136055	MC 30 with printer	30	550 x 640 x 1040	85	3 x 400V, 50/60Hz, 16A, no plug	£35393.00

Accessories

Code	Description	For use with	Price
I136030	Stainless steel cuvette for medium sterilisation	MEDIACLAVE 10	£874.00
I136070	Autoclave kit incl. autoclave cuvette, stainless steel with grid insert and flexible temperature probe	MEDIACLAVE 10	£1254.00
I132130	Magnetic stirrer bar	MEDIACLAVE 10	£211.00
I136034	Decanting tubing incl. silicone tubing, stainless steel nozzle and securing nut	MEDIACLAVE 10	£63.00
I136035	Fitting for dispensing tubing incl. stainless dispense port fitting and spring	MEDIACLAVE 10	£139.00
I136040	Printer splash guard	MEDIACLAVE 10	£48.50
I136247	Injection lid incl. cap, punched disc and septum membrane (silicone/PTFE)	MEDIACLAVE 10	£45.60
I136049	Tubing connector for adding port	MEDIACLAVE 10	£136.00
I136064	Pressure dispensing kit incl. pinch valve box, foot switch, silicone tubing and stainless steel dispensing tube	MEDIACLAVE 10	£810.00
I143200	Foot switch with connecting cable	MEDIACLAVE 10	£67.00
I136060	Stainless steel cuvette incl. stainless steel guide tube for temperature probe	MEDIACLAVE 30	£2132.00
I136075	Magnetic stirrer bar with paddle	MEDIACLAVE 30	£425.00
I136061	Decanting tubing incl. rigid stainless steel tube (length 613 mm), silicone tubing and securing nut	MEDIACLAVE 30	£335.00
I136565	Volume measuring stick for cuvette	MEDIACLAVE 30	£91.00
I136065	Double temperature probe, flexible (PT1000 for MEDIACLAVE, PT100 for external)	MEDIACLAVE 10 & 30	£545.00