

INTEGRA Biosciences MEDIACLAVE Media Preparators Continued

Consumables

Code	Description	Pack	Price
I132150	Silicone lid seal for sealing the vessel lid	Each	£96.00
I136038	Paper rolls for integrated printer	10	£32.50
I136901	Ink ribbon for integrated printer	Each	£30.40
I136047	Septum membrane for injection lid, silicone/PFTE, self-resealing	10	£61.00
I171036	Silicone tubing for pressure dispensing, length 25m (bulk roll), inner Ø 6mm, autoclavable	Each	£157.00



INTEGRA Biosciences MEDIAJET Petri Dish Filler

The unique flexibility to fill Petri dishes of various sizes or with two compartments as well as test tubes of various diameters and length. Truly reliable, walk away operation thanks to a mechanical dish guidance system. Operational downtime and media loss due to 'dish jams' resulting from variations in the disposable plastic Petri dishes are eliminated. In only five minutes, the MEDIAJET vario can be converted from filling 90mm dishes to fill 60 or 35mm dishes. Furthermore, the MEDIAJET includes a shaker mode for facilitating time consuming pour plate applications. The unique TUBEFILLER option is the perfect expansion to the functionality of MEDIAJET, as it allows conversion into a test tube filler in just a minute.

- Fast Petri dish processing for filling up to 900 plates per hour
- Easy conversion at the user level allows the processing of various dish sizes and dish heights
- Biplate option for filling two-compartment Petri dishes with different media
- Shaker mode for pour plate applications
- TUBEFILLER option for the automated filling of up to 1000 test tubes in one hour

Code	Description	Price
I103005	MEDIAJET standard dish filler, base unit only for Ø 90mm Petri dishes, requires rotor and carousel	£25058.00
I103006	MEDIAJET standard dish filler with cooling option, base unit only for 90mm Petri dishes, requires rotor and carousel	£25994.00
I113000	MEDIAJET vario, base unit only for Ø 35, 60 or 90mm Petri dishes, requires rotor, carousel and conversion set	£25898.00
I113001	MEDIAJET vario with cooling, base unit only, for Ø 35, 60 or 90mm Petri dishes, requires rotor, carousel and conversion set	£27061.00
I113002	MEDIAJET vario with biplate & cooling option, base unit only, for Ø 35, 60 or 90mm Petri dishes, requires rotor, carousel and conversion set	£30802.00
I103035	Tubing set for Turbo mode, incl. filling nozzle for Ø 90mm Petri dishes and tubing for double head configuration (6mm ID) - length 2.0m	£325.00
I103036	Turbo mode kit for accelerating the filling rate of Ø 90mm Petri dishes, including extension pump head for double pump head configuration and tubing kit for Turbo mode (Part no. I103035)	£954.00

Carousels and Rotors

An optional inkjet printer is available, please enquire.

Code	Description	For use with	Price
I103020	Carousel, capacity of 360	90mm Ø Petri dishes	£1910.00
I103021	Carousel, capacity of 540	90mm Ø Petri dishes	£3649.00
I113022	Carousel, capacity of 360 (for MJ vario only)	90mm Ø Petri dishes	£1910.00
I113020	Carousel, capacity of 360 (for MJ vario only)	60mm Ø Petri dishes	£2420.00
I113021	Carousel, capacity of 360 (for MJ vario only)	35mm Ø Petri dishes	£3937.00
I103271	Rotor	90mm Ø Petri dishes	£943.00
I103272	Rotor	Small 90mm Ø Petri dishes	£1195.00
I113460	Rotor (for MJ vario only)	Large 90mm Ø Petri dishes	£1195.00
I113271	Rotor (for MJ vario only)	60mm Ø Petri dishes	£1290.00
I113272	Rotor (for MJ vario only)	35mm Ø Petri dishes	£1545.00
I113806	Rotor (for MJ vario with biplate option only)	90mm Ø biplate Petri dishes	£1230.00