

## INTEGRA Biosciences MEDIAJET Petri Dish Filler Continued

**INTEGRA**

### Conversion Sets for MEDIAJET Vario

Code	Description	Comprising	Price
I113051	Conversion set for Ø 90mm petri dishes	Adaption insert, head plates for feeder and stacker piston and 1 tubing set	£1734.00
I113055	Conversion set for large (90L) Ø 90mm petri dishes	Adaption insert, head plates for feeder and stacker piston and 1 tubing set	£1863.00
I113052	Conversion set for Ø 60mm petri dishes	Adaption insert, head plates for feeder and stacker piston, gripper left and 1 tubing set	£2600.00
I113053	Conversion set for Ø 35mm petri dishes	Adaption insert, head plates for feeder and stacker piston, gripper left and gripper right and 1 tubing set	£3236.00

### Accessories

Code	Description	Price
I103030	Tubing set for Ø 90mm Petri dishes, including 1.5m silicone tubing (6 x 9mm) and filling nozzle for Ø 90mm Petri dishes	£253.00
I113030	Tubing set for Ø 60mm and Ø 35mm Petri dishes, including 1.5m silicone tubing (6 x 9mm) and filling nozzle for Ø 60mm and Ø 35mm Petri dishes	£268.00
I103032	Filling nozzle for Ø 90mm Petri dishes	£206.00
I113032	Filling nozzle for 60mm, Ø 35mm Petri dishes (for MEDIAJET vario only)	£217.00
I103040	Tubing set for the preparation of blood agar (suitable to add 3 to 10% blood in connection with the INTEGRA DOSE IT)	£357.00
I103047	Electrical interface cable for combination with DOSE-IT peristaltic pump and MEDIAJET	£172.00
I179147	Suction needle stainless steel	£51.10
I103705	UV lamp TUV11W to reduce the risk of contamination during the filling process	£61.00
I143200	Footswitch	£67.00

### Optional Printer Accessories

Code	Description	Price
I103080	Connecting kit for Imaje 9020 inkjet printer including the fixing device and the interface cable for Ø 90mm and Ø 60mm Petri dishes	£813.00
I103085	Table stand for Imaje 9020 inkjet printer	£1168.00
I113840	Printer connection kit for LINX 4900 inkjet printer for imprinting Petri dishes on the side wall, complete with fixing device and interface cable, without printer	£1278.00
I113841	Cart for MEDIAJET, with compartment for LINX 4900 inkjet printer and opening for connecting the printer head to MEDIAJET	£1664.00
I103091	Printer connection kit for MARKEM-IMAJE 9232 inkjet printer including the fixing device and the interface cable for Ø 90 and Ø 60 mm Petri dishes, without printer	£460.00
I113060	Printer connection kit for Domino A320i inkjet printer including the fixing device and the interface cable for Ø 90 and Ø 60 mm Petri dishes, without printer	£570.00

### Tube-filler Options

Tube-filler includes filling arm, support and a 2.5m (3mm inner diameter) silicone tubing set with dispensing/aspiration and tube collar. Test tube racks, manufactured in polypropylene are autoclavable to 121°C.

**Please note that a centring plate is required to provide accurate positioning of 16, 20 or 25mm test tubes.**

Code	Description	Price
I103010	Tube-filler	£3224.00
I103074	Test tube racks for Ø 30mm test tubes, 3 x 7 capacity	£15.54
I103065	Centring plate for Ø 16mm test tube rack	£234.00
I103066	Centring plate for Ø 20mm test tube rack	£234.00
I103067	Centring plate for Ø 25mm test tube rack	£234.00
I171052	Aspiration/dispensing tube 10cm length, 2.0mm Ø, stainless steel, one end dented	£25.20
I171053	Aspiration/dispensing tube 10cm length, 3.0mm Ø, stainless steel, one end dented	£25.20
I171054	Aspiration/dispensing tube 10cm length, 4.0mm Ø, stainless steel, one end dented	£25.20
I171071	Tube collar as weight for aspiration tubes with 1-3mm inner Ø	£34.80
I171074	Tube collar as weight for aspiration tubes with 4-6mm inner Ø	£28.30