

## OHAUS MB90 and MB120 Moisture Analysers Continued



### Accessories

Code	Alt Ref	Description	Pack	Price
BAL5536	30246947	Pan holder	Each	£41.40
BAL5538	30284478	In-use cover	Each	£35.40
BAL5540	30284477	Scoop	Each	£28.50
BAL6240	11113857	Temperature calibration kit	Each	£309.00
BAL5152	80850043	Security device	Each	£157.50
BAL4237	30585411	Aluminum pans	50	£13.20
BAL6236	80850087	Glass fiber filters	200	£50.30
BAL5542	80252477	Sample cage	Each	£41.00
BAL5544	80252478	Reusable pans, 7mm	3	£86.30
BAL5546	80252479	Reusable pans, 14mm	3	£88.60
BAL1021	30064202	Printer SF40A	Each	£466.00

## OHAUS MB23 Moisture Analyser



The MB23 provides precise moisture content determination at an unbeatable value. It uses glassless infrared heating, which increases efficiency for faster results. With a large backlit LCD display, standard RS232 port, 110g capacity and a readability of 0.01g/0.1% the MB23 offers moisture analysis for any budget.

- Glassless infrared heating increases efficiency for faster results
- User friendly touch button design for easy operation
- Compact design suits any workspace



Code	Alt Ref	Model	Capacity, g	Temp Range, °C	Readability, mg	Pan Dims, mm	Price
BAL5502	80252470	MB23	110	50 - 160 (increments of 5)	10/0.1%	90	£1150.00

### Accessories

Code	Alt Ref	Description	Pack	Price
BAL4237	30585411	Aluminum pans	50	£13.20
BAL6234	80850088	Reusable sample pans	3	£86.30
BAL6236	80850087	Glass fiber filters	200	£50.30
BAL6240	11113857	Temperature calibration kit	Each	£309.00

## Can't find what you're looking for?

Let us know! Our team are available to help. Contact your local Territory Sales Manager to discuss your requirements.



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

