



Kern F1 Weights

KERN

Compact, individual weights, made of polished stainless steel. Includes container manufactured from lined plastic and DKD certificate.

Code	Alt Ref	Weight, g	Tolerance, ± mg	Price
BAL1776	326-014-600	1	0.10	£49.00
BAL1778	326-024-600	2	0.12	£58.00
BAL1780	326-034-600	5	0.16	£55.00
BAL1782	326-044-600	10	0.20	£58.00
BAL1784	326-054-600	20	0.25	£59.00
BAL1786	326-064-600	50	0.30	£71.00
BAL1788	326-074-600	100	0.50	£68.00
BAL1790	326-084-600	200	1	£119.00
BAL1792	326-094-600	500	2.50	£121.00
BAL1794	326-114-600	1000	5	£154.00
BAL1796	326-124-600	2000	10	£260.00



Kern F1 Weight Sets

KERN

ECO-shape weights supplied in a plastic case. Milligram weights manufactured from aluminium/German silver; individual weights made from polished stainless steel. Milligram weight sets 1 - 500mg are supplied in a separate plastic case. Available with or without DAKS calibration certificate.

Code	Alt Ref	Nominal Weight, g	Calibration Certificate	Price
BAL3152	325-024	0.001 - 50	No	£729.00
BAL3154	325-024+962-401	0.001 - 50	Yes	£661.00
BAL3156	325-034	0.001 - 100	No	£782.00
BAL3158	325-034+962-402	0.001 - 100	Yes	£679.00
BAL3182	325-044	0.001 - 200	No	£902.00
BAL3184	325-044+962-403	0.001 - 200	Yes	£802.00
BAL3186	325-054	0.001 - 500	No	£1022.00
BAL3188	325-054+962-404	0.001 - 500	Yes	£880.00
BAL3246	326-024	1 - 50	No	£489.00
BAL3248	326-024 + 962-415	1 - 50	Yes	£352.00
BAL3282	326-034	1 - 100	No	£542.00
BAL3284	326-034+962-416	1 - 100	Yes	£385.00
BAL3286	326-044	1 - 200	No	£684.00
BAL3288	326-044 + 962-417	1 - 200	Yes	£527.00
BAL3306	326-054	1 - 500	No	£782.00
BAL3308	326-054 + 962-418	1 - 500	Yes	£596.00

View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

scientificlabs.co.uk/literature (UK & Northern Ireland)

scientificlabs.ie/literature (Republic of Ireland)

scientificlabs.com/literature (International)

