


Thermo Scientific™ Nunc™ 384 Well Microplates

Thermo
SCIENTIFIC

- Improve signal to noise ratio by minimising background autofluorescence and well to well crosstalk
- Black microplates are recommended for fluorescence measurements, with minimum backscattered light and background fluorescence
- White microplates are recommended for luminescence measurements, with maximum signal reflection and minimal crosstalk between wells
- Low flange skirt

Code	Description	Colour	Material	Sterility	Pack	Price
142761	Tissue culture treated, optical bottom	Black	Polymer	Sterile	30	£694.00
142762	Tissue culture treated, optical bottom	White	Polymer	Non sterile	30	£784.00
152029	Poly-D-Lysine coated, optical bottom	Black	Polymer	Non sterile	20	£1311.00
164564	Cell culture treated	Black	Polystyrene	Sterile	30	£408.00
164688	Cell culture treated	Clear	Polystyrene	Sterile	30	£319.00
242764	Non treated	Black	Polymer	Non sterile	30	£340.00
262160	Non treated	Clear	Polystyrene	Non sterile	100	£598.00
262260	Non treated	Black	Polystyrene	Non sterile	100	£799.00
264704	Non treated, shallow well	Clear	Polystyrene	Non sterile	100	£821.00
264705	Non treated, shallow well	Black	Polystyrene	Non sterile	100	£325.00




Hellma Quartz Microplates

Hellma
Analytics

- 96 well
- Synthetic quartz glass
- 2mm base thickness.
- Resistant to aggressive fluids
- Optically suitable for microscopy
- Suitable for UV and fluorescence applications
- 2mm base thickness

Code	Alt Ref	Well Volume, µL	Dims, w x d x h, mm	Price
CEL1732	730-009-44	30	127 x 85.5 x 14.5	£3672.00



SLS Select ClearPeel Heat Seal Microplate Film

SLS
SELECT

- Peelable optically clear heat sealing film
- Suitable for qPCR and optical applications.
- Optically clear laminate film forming a peelable seal to polypropylene, polyethylene, polystyrene, polycarbonate and cyclic olefin copolymer (COC) plates
- Seal integrity range: -80 to 80°C (to 110°C with a pressurised heated PCR lid)
- Moderate solvent resistance
- Suitable for short term compound storage at room temperature

Code	Dims, w x l, mm	Pack	Price
SLS9101	78 x 125	100	£43.10
SLS9100	78 x 500000	Roll	£425.00