

Cerillo Stratus Mini Plate Readers Continued

Adapter

The adapter securely attaches the Stratus to a wide variety of small orbital shakers and shaking incubators.


Code	Description	Price
MIC6036	Adapter	£375.00

Tecan HydroFlex Plate Washers



Provides excellent automated microplate strip washing and vacuum filtration performance for 96 well microplate formats. This modular and upgradeable platform is ideal for a wide range of cell, enzyme and DNA based applications in academia, biotech, pharma and clinical diagnostics, reflecting over 30 years of Tecan expertise in advanced liquid handling. The multifunctional HydroFlex microplate washer offers far greater flexibility than a standard plate microplate washer, and features user-interchangeable plate carriers for easy switching between washing and filtration applications to allow:

- ELISA washing for improved consistency and reliability
- Gentle washing of cell based assays with adherent or weakly adherent cells
- Microplate array washing for protein and antibody arrays
- Vacuum filtration to waste for PCR cleanup after DNA amplification, etc
- Magnetic bead washing using MBS carrier
- Compact and modular design, suitable for a range of applications and laboratories
- Ready for automated washing of ELISAs, cells and protein arrays
- Automated system for vacuum filtration to waste
- Magnetic plate carrier for magnetic bead washing
- Advanced online process controls for safety and reliability
- Bubble sensor unit that ensures reliable buffer dispensing
- Multipoint aspiration for flat bottom plates to achieve minimal residual volumes
- Individual software control of speed settings and wash head positions
- Suitable for integration onto Freedom EVO and Freedom EVOLyzer workstations
- Excellent reliability and low service costs
- Level detection option available, please contact customer services for more information



1
YEAR

Code	Alt Ref	Channel	No. of Inlets	Price
MIC9288	30082285/30023426	8	1	£4819.00
MIC9290	30082285/30023425	8	2	£5399.00
MIC9292	30082285/30023424	8	4	£5541.00
MIC9294	30022337/30023426	16	1	£5297.00
MIC9296	30022337/30023425	16	2	£5604.00
MIC9298	30022337/30023424	16	4	£5763.00

Accessories

Code	Alt Ref	Description	Price
MIC9300	30023217	MBS plate carrier with 96 magnets (order with instrument only)	£1159.00
MIC9302	30034052	Smart-MBS plate carrier with 108 magnets	£1811.00
MIC9304	30041515	Smart 2-MBS plate carrier with 108 magnets	£2503.00
MIC9306	30032418	Low level vacuum filtration incl. plate carrier with soft-gasket (order with instrument only)	£4137.00
MIC9308	30024633	High level vacuum filtration incl. plate carrier with hard-gasket (order with instrument only)	£7241.00