

SLS Select Weighing Boats

Anti static polystyrene weighing boats.



White

Code	Capacity, mL	Dims, w x d, mm	Shape	Pack	Price
BAL1818	5	35 x 55	Diamond	1000	£39.60
BAL1820SP	7	44 x 44	Square	50	£2.12
BAL1820	7	44 x 44	Square	1000	£30.50
BAL1826	30	65 x 85	Diamond	500	£52.00
BAL1822SP	100	80 x 80	Square	50	£3.20
BAL1822	100	80 x 80	Square	1000	£43.80
BAL1828	100	100 x 125	Diamond	250	£36.10
BAL1824SP	250	134 x 134	Square	50	£5.20
BAL1824	250	134 x 134	Square	1000	£96.00



Black

Code	Capacity, mL	Dims, w x d, mm	Shape	Pack	Price
BAL1819	5	35 x 55	Diamond	1000	£45.70
BAL1821	7	44 x 44	Square	1000	£51.00
BAL1827	30	65 x 85	Diamond	1000	£35.00
BAL1823	100	80 x 80	Square	1000	£46.40
BAL1829	100	100 x 125	Diamond	500	£43.00
BAL1825	250	140 x 140	Square	250	£115.00



Aluminium Dishes

Inexpensive, pure aluminium dishes are suitable for the processing/testing of both wet and dry samples. They enable wet solids testing from mixtures and emulsions by the loss of weight method without the need of time-consuming scrubbing of other more expensive dishes.

- Use once and discard
- Unaffected by thermal shock
- Suitable for temperatures up to 450°C
- Can be used in pairs for covered samples
- With flange

Code	Dims, h x Ø, mm	Pack	Price
DIS5000	5 x 42	500	£132.00
DIS5004	10 x 55	250	£119.00



Heathrow Scientific Aluminium Moisture Dishes



Designed for moisture determination balances. Slightly tapered walls for easy stacking.

Code	Alt Ref	Dims, h x Ø, mm	Pack	Price
MOI2101	HS14522	8 x 101	50	£42.90