

Heathrow Scientific Aluminium Weighing Boats



Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling and a flanged, crimped wall for rigidity. Finger tab handle can be etched for identification prior to weighing.

Supplied in packs of 100

Code	Alt Ref	Description	Dims, h x Ø, mm	Price
MOI4004	HS14521A	Small	12 x 43	£25.00
MOI4002	HS14521B	Medium	14 x 57	£27.20
MOI4000	HS14521C	Large	16 x 70	£32.50



PYREX® Glass Weighing Boats



Ideal for weighing small quantities of powder. Hollow stem aids the introduction of powder into receiving vessel, and a flat base ensures stability on the balance pan.

Code	Alt Ref	Capacity, mL	Length, mm	Price
BAL1870	2115/02M	3	70	£15.80
BAL1872	2115/04M	6	85	£17.38
BAL1874	2115/06M	10	100	£20.93



Glass Weighing Bottles

Weighing bottles with ground hollow blown and totally enclosed stopper. Made of borosilicate glass 3.3, DIN 12605.

Code	Description	Dims, h x Ø, mm	Price
BAL1852	Squat	30 x 60	£42.00
BAL1862	Tall	60 x 30	£16.00



Whatman Weighing Paper



- Grade B-2
- Dissolves residue free in acid solution, which is critical in applications such as Kjeldahl analysis of nitrogen in soil
- Thin (0.04mm) and lightweight (43g/m²)
- Transparent and smooth parchment surface
- Simple weighing out and sample transfer
- Reliable for quantitative transfer
- Low nitrogen paper that is residue free for accurate chemical analysis

Supplied in packs of 500 sheets

Code	Dims, w x l, mm	Price
Z134112-500EA	76.2 x 76.2	£31.00
BAL1882	101.6 x 101.6	£29.00
BAL1884	152.4 x 152.4	£30.00

