



Watch Glasses

- Soda lime glass
- Ground edges and round bottom

Supplied in packs of 10

Code	Ø, mm	Price
WAT1004	50	£10.60
WAT1006	60	£11.50
WAT1008	70	£11.80
WAT1010	80	£11.80
WAT1012	90	£13.20
WAT10140	100	£11.80
WAT10160	130	£17.20
WAT10180	150	£23.70
WAT10200	200	£42.40




PYREX® Crystallising Basins

PYREX®

- Flat bottomed with straight sides and no spout
- Complies to DIN 12337
- Manufactured from PYREX® borosilicate glass
- Suitable for use up to 400°C
- Chemically and thermally resistant

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BAS1518	1470/03M	15	40 x 25	10	£90.00
BAS1516	1470/01M	60	58 x 35	10	£72.00
BAS15200	1470/02D	100	70 x 40	10	£79.00
BAS15220	1470/04D	150	80 x 45	10	£97.00
BAS15240	1470/08D	300	95 x 55	10	£99.00
BAS15260	1470/10D	500	115 x 65	10	£172.00
BAS15280	1470/12D	900	140 x 75	10	£234.00
BAS1530	1470/16D	2000	190 x 90	Each	£48.80



PYREX® Watch Glasses

PYREX®

- Manufactured from PYREX® borosilicate glass
- High resistance to thermal shock and chemical attack
- For evaporation work and covering heated beakers
- Fire polished edges to reduce chipping and breakage
- Complies with DIN 12341

Supplied in packs of 10

Code	Alt Ref	Ø, mm	Price
WAT10500	1060/04D	40	£117.00
WAT10520	1060/06D	50	£117.00
WAT10540	1060/09D	60	£129.00
WAT10560	1060/10D	80	£134.00
WAT10580	1060/14D	100	£151.00
WAT10600	1060/18D	125	£171.00
WAT10620	1060/20D	150	£205.00
WAT10610	1060/21D	200	£313.00