



### Motic SMZ-161 Series Stereo Microscopes

**Motic®**

Convincing optical performance, combined with a small footprint and robust mechanics, makes the SMZ-161 the ideal stereo microscope for teaching in laboratories. Greenough zoom optical system, 45°, 360° rotating inclined head with widefield eyepieces WF10X/20mm with diopter adjustment on both eyepiece tubes (50 – 75mm), zoom system ratio 6:1, standard working distance 110mm, 3W LED with intensity control for incident and transmitted illumination, black/white and frosted stage plates and dust cover.

Code	Alt Ref	Description	Price
<b>MIC1588</b>	1100201300062	SMZ-161-binocular LED	£835.00
<b>MIC1590</b>	1100201300082	SMZ-161-trinocular LED	£897.00



### Motic SMZ-171 Series Stereo Microscopes

**Motic®**

A versatile stereo microscope platform for biological and medical preparations to industrial quality control and inspection. Greenough zoom optical system, 45° inclined binocular or trinocular head with widefield adjustable eyepieces N-WF10X/23mm with diopter adjustment, ±5 diopter on both eyepiece tubes (52 – 75mm), zoom system ratio 6.7:1, standard working distance 110mm, transmitted light with tilting reflector (oblique illumination). Black/white and frosted glass stage plates and dust cover.

Code	Alt Ref	Description	Price
<b>MIC5444</b>	1100200600754	SMZ-171-binocular LED	£1400.00
<b>MIC5445</b>	1100200600766	SMZ-171-trinocular LED	£1484.00
<b>MIC5446</b>	1100200600758	SMZ-171-binocular plain stand	£1092.00
<b>MIC5447</b>	1100200600762	SMZ-171-trinocular plain stand	£1174.00

#### Accessories and Spares

Code	Alt Ref	Description	Price
<b>MIC5448</b>	1101001703821	Objective 0.5x ESD (WD=191.8mm)	£115.00
<b>MIC5449</b>	1101001703661	Objective 0.63x ESD (WD=142.7mm)	£115.00
<b>MIC5450</b>	1101001703811	Objective 0.75x ESD (WD=128.6mm)	£115.00
<b>MIC5451</b>	1101001703671	Objective 1.5x ESD (WD=56.3mm)	£115.00
<b>MIC5452</b>	1101002402721	LED ring illumination 60T dimmable and segmentable	£202.00
<b>MIC2806</b>	1101002000011	MLC-150 fibre optic illuminator	£495.00
<b>MIC2804</b>	1101002000142	Gooseneck bifurcated light guide	£212.00
<b>MIC5455</b>	1101002000072	Blue filter for MLC-150	£41.00

## Authorised distributor



We offer a broad portfolio of products from key suppliers such as:

**Reagecon**  
Delivering the Correct Result

**Honeywell**

**MERCK**

**Sigma-Aldrich®**  
Lab & Production Materials

**Supelco®**  
Analytical Products

**spex®**  
an anlyia scientific company