

## Olympus SZ Series Stereo Microscopes

**OLYMPUS**

The SZ series of stereo microscopes offer a variety of options to cover life science applications. Available with or without illumination and in binocular or trinocular setup, the SZ series Greenough optical system offer a small angle convergence of the image forming path for high levels of flatness and depth of focus. Optical coatings ensure high colour fidelity whilst antistatic materials and coatings protect specimens from anti static discharge. LED illumination versions include the option of LED transmitted light with mirror contrast or and LED stand with simultaneous transmitted and reflected light.

- ComfortView eyepieces offer pupil aberration control
- Optical coating ensure high quality colour and protect specimens from anti-static discharge
- Options include no illumination or LED illumination via mirror contrast or an LED stand
- Available in binocular or trinocular formats
- Olympus microscopes are supplied with 12 month on site Ultimate Service Care Plan – see start of microscope listing for details



### SZ Series with Basic Stand and No Illumination

- Binocular or trinocular microscope body
- 0.5x C-mount camera adapter on trinocular version
- Standard stand
- Eyepiece
- Stage plate

Code	Alt Ref	Description	Zoom Ratio	Magnification	Price
<b>MIC1086</b>	K23008479	SZ51 binocular with basic stand	5:1	0.8x to 4x	<b>£1842.00</b>
<b>MIC1092</b>	K23008482	SZ61 trinocular with basic stand	6.7:1	0.67x to 4.5x	<b>£4589.00</b>

### SZ Series with Basic Stand and LED Transmitted Light (with Mirror Contrast Adjustment)

- Binocular or trinocular microscope body
- 0.5x C-mount camera adapter on trinocular version
- Standard stand
- Eyepiece
- Stage glass
- Light guide adapter
- Focusing lens
- LED light source with tiltable mirror providing direct and oblique illumination
- Flexible light guide

Code	Alt Ref	Description	Zoom Ratio	Magnification	Price
<b>MIC1088</b>	K23008480	SZ51 binocular with basic stand - LED transmitted light	5:1	0.8x to 4.5x	<b>£2762.00</b>
<b>MIC1094</b>	K23008483	SZ61 trinocular with basic stand - LED transmitted light	6.7:1	0.67x to 4.5x	<b>£5451.00</b>

### SZ Series with LED Base for Transmitted (Fixed Angle) and Reflected (Variable Angle) Light

- Binocular or trinocular microscope body
- 0.5x C-mount camera adapter on trinocular version
- Eyepiece
- LED base with transmitted and reflected light LED illumination for brightfield observation

Code	Alt Ref	Description	Zoom Ratio	Magnification	Price
<b>MIC1980</b>	K23008952	SZ51 binocular with LED base	5:1	0.8x to 4x	<b>£2938.00</b>
<b>MIC1096</b>	K23008484	SZ61 trinocular with LED base	6.7:1	0.67x to 4x	<b>£5759.00</b>