

SAMCO Forceps



- Stainless steel
- Spatulate tips for lifting microscopical sections, slides or coverslips



Code	Alt Ref	Price
INS4420	L325/01	£5.90

Heathrow Scientific Tall Slide Box



- For storage and transport of larger, non standard microscope slides (3 x 2")
- Made from rugged ABS materials for safe storage and transport
- Grooved slots hold individual slides in place
- Soft cork lining insures that slides will not slip
- Identification card, lining and grooved slots are numbered
- Stackable with other tall boxes and/or 100 place premium slide boxes



Code	Alt Ref	Capacity	Colour	Dims, l x w x h, mm	Price
MIC4080	HS15979A	100 slides	Blue	208 x 175 x 60	£32.10

Heathrow Scientific Economy Slide Boxes



- Hold standard slides (75 x 25mm and 3 x 1")
- Constructed of durable polypropylene
- Hinged lids are held securely closed by latch lock
- Stacking feet maximise stability
- Slide identification card
- No liner



Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm	Price
MIC3580	HS15990A	25 slides	Blue	141 x 36 x 92	£11.21
MIC3582	HS15990B	25 slides	Red	141 x 36 x 92	£11.21
MIC3584	HS15990C	25 slides	White	141 x 36 x 92	£11.80
MIC3586	HS15991A	100 slides	Blue	169 x 37 x 210	£15.00
MIC3588	HS15991B	100 slides	Red	169 x 37 x 210	£15.00
MIC3590	HS15991C	100 slides	White	169 x 37 x 210	£14.10

Heathrow Scientific 25 Place Premium Slide Boxes



- Hold standard slides (75 x 25mm and 3 x 1")
- Moulded, rigid boxes
- Securely hold slides in grooved slots
- Durable ABS plastic with cork lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet
- Easily access slides with friction fit cover



Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm	Price
MIC3562	HS15989A	25 slides	Blue	88 x 141 x 35	£16.10
MIC3564	HS15989B	25 slides	Green	88 x 141 x 35	£16.10
MIC3566	HS15989C	25 slides	Red	88 x 141 x 35	£16.10
MIC3570	HS15989E	25 slides	White	88 x 141 x 35	£16.10