



Heathrow Scientific 50 Place Premium Slide Boxes



- Hold standard slides (75 x 25mm and 3 x 1")
- Moulded, rigid boxes
- Securely hold slides in grooved slots
- Durable ABS plastic with cork lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet
- Thumb latch lock mechanism

Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm	Price
MIC4214	HS15996A	50 slides	Blue	86 x 209 x 35	£17.91
MIC4216	HS15996B	50 slides	Red	86 x 209 x 35	£16.44
MIC4004	HS15996C	50 slides	White	86 x 209 x 35	£9.70



Heathrow Scientific 100 Place Premium Slide Boxes



- Hold standard slides (75 x 25mm and 3 x 1")
- Moulded, rigid boxes
- Securely hold slides in grooved slots
- Durable ABS plastic with cork lining
- Both the lining and the slots are numbered to correspond to the slide inventory sheet
- Stay organised with a moulded labelling area with "Project" and "Date"
- Closes securely with a rust resistant nickel plated clasp and pin hinges

Code	Alt Ref	Capacity	Colour	Dims, w x d x h, mm	Price
MIC3552	HS15994B	100 slides	Green	175 x 208 x 34	£18.20
MIC3560	HS15994F	100 slides	Grey	175 x 208 x 34	£22.50
MIC4012	HS15994G	100 slides	Black	175 x 208 x 34	£20.40



Heathrow Scientific Microscope Slide Mailers



- Polypropylene containers which hold five standard slides (25 x 75 mm and 1 x 3")
- Ideal for shipping, storing, or staining slides
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe
- Disposable

Supplied in packs of 25

Code	Alt Ref	Style	Colour	Dims, l x w x h, mm	Price
MIC4040	HS15982	Side opening	Natural	81 x 19 x 31	£38.70
MIC4042	HS15983G	End opening	Red, yellow, blue, purple, orange	43 x 24 x 88	£41.40



Cole-Parmer® Microscope Slide Warmer



- Facilitates tissue mounting and slide staining procedures
- Dehydrate and set adhesives or perform non adhesive section mountings on warming deck
- Front panel has temperature adjustment dial and heater/power on lights
- Drying bars included for easier slide drying

Code	Alt Ref	Model	Capacity	Temp Range, °C	Dims, w x d x h, mm	Weight, kg	Price
MIC4302	48573-06	MH6616	50 slides	Ambient to 100	394 x 181 x 89	1.8	£451.00