



## Thermo Scientific™ Vortex Mixer

**Thermo**  
SCIENTIFIC

Designed for accurate speed control for continuous mixing or pulsing of samples.

- Digital display of speed and time
- Continuous and touch modes
- Precise speed control:  $\pm 5\%$  RPM
- Speed accuracy  $\pm 5\%$  <1500RPM,  $\pm 10\%$  >1500RPM
- Rapid acceleration to maximum RPM, touch: 1500RPM in 3 seconds, continuous: 3000RPM in 6 seconds
- Continuous or countdown timer (up to 99 minutes)
- Maximum load of 0.5kg at 800RPM
- Includes rubber platform

Code	Alt Ref	Speed Range, RPM	Orbit, mm	Dims, w x d x h, mm	Weight, kg	Price
MIX1312	88882010	200 - 3000	4	151 x 205 x 131	6.6	£420.00

### Accessories

Code	Alt Ref	Description	Mode	Max Speed, RPM	Price
MIX1364	88882120	Rubber platform	Touch/continuous	3000	£18.00
MIX1366	88882121	Flat bottom vessel platform	Touch/continuous	3000	£18.00
MIX1368	88882122	Microplate tray	Continuous	1200	£18.00
MIX1370	88882123	Tube holder 10 x 19mm	Continuous	2000	£25.00
MIX1372	88882124	Tube holder 12 x 13mm	Continuous	1800	£25.00
MIX1374	88882125	Tube holder 15 x 12mm	Continuous	1000	£25.00
MIX1376	88882126	Tube holder 20 x 7mm	Continuous	1000	£25.00
MIX1378	88882127	Tube holder 26 x 4mm	Continuous	900	£25.00



## Benchmark Vornado™ Mini Vortex Mixers

**Benchmark**  
Scientific

Powerful is the best word to describe this mini vortexer. With a 4mm orbit and fixed speed of 2800RPM, the Vornado™ instantly vortexes even the largest samples, including nearly full 50mL tubes. Despite its powerful motor, the Vornado™ has a footprint of less than 4 x 4 inches, allowing it to fit on even the most crowded bench. Insta-Touch operation means that the motor is automatically activated with slight pressure on the centre cup head. Operation is smooth, quiet and efficient. When pressure is removed, operation immediately ceases. The Vornado™'s unique head design prevents liquids from entering the housing, prolonging motor life.

- Powerful vortexing for tubes up to 50mL
- Inst-Touch™ pressure activated operation
- "Mini sized", less than 4" wide
- Speed: 2800RPM
- Orbit: 4mm
- Dims, w x d x h: 94 x 99 x 66mm
- Weight: 0.4kg

Code	Alt Ref	Colour	Price
MIX7000	BV101-E	Grey	£211.00
MIX7002	BV101-B-E	Blue	£211.00
MIX7004	BV101-G-E	Green	£211.00
MIX7006	BV101-P-E	Purple	£211.00
MIX7008	BV101-R-E	Red	£211.00