

## Eppendorf Innova® Stackable Incubated/Refrigerated Shakers

eppendorf

Heavy duty triple eccentric drive providing uniform motion and superb reliability. All units have audible and visual alarms, operate at 25 to 400RPM, timer (0-99.9 hrs) with continuous operation or automatic shut off at the end of the study.

- Temperature range: ambient +5 to 80°C
- Refrigerated versions (R) temperature range: ambient -20 to 80°C (minimum set point is 4°C)
- 44/44R shakers can be stacked up to 3 high and Innova® 42/42R shakers can be stacked up to 2 high
- 42/42R models have been designed to use with a maximum flask size of 4L, whereas 44/44R models can accommodate 5L flasks
- Platforms are sold separately
- **Note:** see page 680 for platforms – Innova® 2100 platforms are the same ones used in the 42/42R and 44/44R universal platform is listed in the accessories below; please contact Customer Services for details of other platforms



Code	Alt Ref	Model	Dims, w x d x h, mm	Price
SHA4135	M1335-0002	42	620 x 754 x 825	£10373.00
SHA4134	M1335-0012	42	620 x 754 x 825	£10373.00
SHA4139	M1335-0006	42R	620 x 754 x 825	£11777.00
SHA4140	M1335-0016	42R	620 x 754 x 825	£11777.00
SHA1100	M1282-0002	44	1350 x 870 x 735	£14586.00
SHA1102	M1282-0012	44	1350 x 870 x 735	£14586.00
SHA1104	M1282-0006	44R	1350 x 870 x 735	£17945.00
SHA1106	M1282-0016	44R	1350 x 870 x 735	£17945.00

### Options/Accessories

Code	Alt Ref	Description	Price
SHA1109	M1335-0800	Stacking kit for 2 x Innova® 42/42R	£361.00
SHA1101	M1335-0810	Stacking kit to stack Innova® 42/42R units on Innova® 4200/4230	£517.00
SHA1110	M1282-0600	Short base to raise shakers by 102mm, for use with up to 3 stacked units	£659.00
SHA1112	M1282-0800	Medium base to raise shakers by 330mm, for use with up to 2 stacked units	£1208.00
SHA1114	M1282-0700	Tall base to raise shakers by 432mm, for use with a single unit	£1309.00
SHA1116	M1282-0500	Stacking kit for 2 x Innova® 44/44R	£362.00
SHA1128	M1282-9921	UV germicidal lights (factory installed)	£1992.00
SHA1122	M1282-9922	Photosynthetic light bank (for Innova® 44R)	£1992.00
SHA1124	M1282-0291	6 port gassing unit for direct gassing into flasks (factory installed)	£737.00
SHA1126	M1282-5011	Humidity monitor (factory installed)	£998.00
SHA4070	M1250-9902	Universal platform for Innova® 42	£820.00
SHA1108	M1282-9904	Universal platform for Innova® 44	£1150.00

### KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty