



Eppendorf Test Tube Racks

eppendorf

Test tube racks with an adjustable angle and tube height support. Manufactured from stainless steel with inserts to help retain tubes and prevent them from spinning in the rack. Finger tight knobs let you quickly adjust the angle, or release the rack from the base plate, for easy transport to the bench.

Small Racks

Code	Alt Ref	Tube Ø, mm	No. of Places	Price
SHA4400	M1289-0001	8 - 11	48	£319.00
SHA4402	M1289-0002	12 - 15	34	£319.00
SHA4406	M1289-0004	18 - 21	18	£319.00
SHA4408	M1289-0005	22 - 26	13	£319.00
SHA4452	M1289-0006	26 - 30	12	£319.00

Medium Racks

Code	Alt Ref	Tube Ø, mm	No. of Places	Price
SHA4410	M1289-0010	8 - 11	60	£319.00
SHA4412	M1289-0020	12 - 15	44	£319.00
SHA4414	M1289-0030	15 - 18	31	£319.00
SHA4416	M1289-0040	18 - 21	23	£319.00
SHA4418	M1289-0050	22 - 26	16	£319.00
SHA4454	M1289-0060	26 - 30	16	£319.00

Large Racks

Code	Alt Ref	Tube Ø, mm	No. of Places	Price
SHA4420	M1289-0100	8 - 11	80	£319.00
SHA4456	M1289-0200	12 - 15	60	£319.00
SHA4458	M1289-0300	15 - 18	42	£319.00
SHA4460	M1289-0400	18 - 21	30	£319.00
SHA4462	M1289-0500	22 - 26	22	£319.00
SHA4464	M1289-0600	26 - 30	20	£319.00



Cole-Parmer® Stuart™ Incubated Shakers

Cole-Parmer
stuart™

These shaking incubators provide an orbital shaking motion to ensure even mixing while promoting bacterial growth. All models feature a timer which can be set from 1 second to 9 days. Once the timer has counted down, an alarm will sound and the shaking action will stop while the incubator will continue to run. The retractable platform allows easy access to the samples at the back of the incubator yet securely locks in place during operation. Built in USB connection and dedicated software allow long term monitoring of the incubator temperature via a PC.

- Shaker and incubator combination makes these ideal for cell culture procedures
- Microprocessor controlled shaking speed ensures accuracy
- Forced air circulation ensures uniform temperature inside the chamber
- Easily set and view temperature and speed on the digital display
- Unique retractable platform offers easy loading and unloading

Code	Alt Ref	Model	Capacity, L	Temp Range, °C	Platform Dims, mm	Overall Dims, w x d x h, mm	Weight, kg	Heat Output, W	Price
SHA3010	79520-02	SI500	51	Ambient +5 to 60	335 x 335	450 x 474 x 522	30	250	£3915.00
SHA3006	79520-16	SI600	115	Ambient +5 to 60	526 x 390	675 x 542 x 642	55	400	£5187.00
SHA3008	79520-22	SI600C	115	Ambient -15 to 60	520 x 390	675 x 562 x 640	60	400	£6202.00