

## Grant JB Academy Unstirred Water Baths

grant<sub>o</sub>

A basic, simple to use, high quality water bath range, consisting of three great value models. A polycarbonate, perforated base tray is included as standard (lid not included), which promotes heat convection, optimal temperature uniformity, and allows any container to be placed directly in the bath.

- Temperature range: ambient +5 to 95°C ( $\pm 0.5^\circ\text{C}$ )
- Unique Set and Forget<sup>®</sup> technology with fast heat up
- Practical front panel lock which disables front panel controls preventing unintentional temperature changes



Code	Alt Ref	Capacity, L	Ext Dims, w x d x h, mm	Int Dims, w x l, mm	Drain Outlet	Price
<b>BAT3310</b>	JBA5	5	335 x 215 x 205	281 x 131	No	£389.00
<b>BAT3312</b>	JBA12	12	360 x 390 x 227	306 x 281	No	£448.00
<b>BAT3314</b>	JBA18	18	335 x 590 x 277	281 x 485	No	£520.00

## Grant Lids for Unstirred Water Baths

grant<sub>o</sub>

Use lids to direct condensation away from immersed vessels, avoid contamination, reduce evaporation and save energy. Recommended when using water baths above 60°C. Flat lids have ring sets of variable hole diameter to accommodate tall vessels.



Code	Description	For use with	Price
<b>BAT3031</b>	Gabled polycarbonate lid	SAP2	£58.00
<b>BAT3036</b>	Gabled polycarbonate lid	JBN5, JBA5, SAP2S, SAP5, SBBAQUA5+	£73.00
<b>BAT3037</b>	Gabled polycarbonate lid	JBN12, JBA12, SAP12, SBBAQUA12+	£89.00
<b>BAT3039</b>	Gabled polycarbonate lid	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA26+	£98.00
<b>BAT3120</b>	Gabled stainless steel lid	JBN5, JBA5, SAP2S, SAP5, SBBAQUA5+	£138.00
<b>BAT3122</b>	Gabled stainless steel lid	JBN12, JBA12, SAP12, SBBAQUA12+	£151.00
<b>BAT3124</b>	Gabled stainless steel lid	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA18+, SBBAQUA26+	£173.00
<b>BAT3126</b>	Gabled stainless steel lid	SAP34	£199.00
<b>BAT3130</b>	Flat stainless steel lid with 2 ring sets	JBN5, JBA5, SAP5, SBBAQUA5+	£167.00
<b>BAT3132</b>	Flat stainless steel lid with 4 ring sets	JBN12, JBA12, SAP12, SBBAQUA12+	£256.00
<b>BAT3134</b>	Flat stainless steel lid with 6 ring sets	JBN18, JBN26, JBA18, SAP18, SAP26, SBBAQUA18+, SBBAQUA26+	£362.00
<b>BAT3136</b>	Flat stainless steel lid with 8 ring sets	SAP34	£418.00