



Simax® Glass Mortars and Pestles

SIMAX®

Manufactured using borosilicate glass PN 70 0080 with even and consistent wall thickness to ensure the longevity. Glazed inside.

Code	Alt Ref	Mortar Ø, mouth x base, mm	Mortar Height, mm	Price
MOR1000	2231/80	80 x 50	60	£24.40
MOR1002	2231/100	100 x 60	75	£32.30
MOR1004	2231/150	150 x 85	110	£54.40



Agate Mortar and Pestle

Highly polished natural stone with working surfaces free from blemishes. Grinding surface unglazed.

Code	Mortar Ø, Ext x Int, mm	Price
MOR4000	50 x 40	£143.00



Porcelain Mortars and Pestles

Unglazed grinding surface.

Code	Mortar Capacity, mL	Mortar Dims, h x Ø, mm	Pestle Dims, l x Ø, mm	Price
MOR3000	30	36 x 56	115 x 24	£15.00
MOR3002	75	45 x 76	115 x 24	£17.20
MOR3004	160	56 x 90	135 x 30	£23.50
MOR3006	220	64 x 105	150 x 36	£27.80
MOR3008	400	71 x 125	150 x 36	£46.50
MOR2010	650	90 x 150	210 x 55	£66.00



Mortars and Pestles, Porcelain

Unglazed grinding surface.

Code	Capacity, mL	Mortar Ø, mm	Price
PM100	100	100	£9.04
PM90	150	90	£10.22
PM130	300	130	£11.52